

## Sample Preparation for X-ray Micro-computed Tomography of Woody Plant Material and Associated Xylem Visualisation Techniques

Nicholas J. B. Brereton\*

Institut de recherche en biologie végétale, Université de Montréal, Montreal, Canada

\*For correspondence: [Nicholas.brereton@umontreal.ca](mailto:Nicholas.brereton@umontreal.ca)

**[Abstract]** Variation in the tissue structure of short rotation coppice (SRC) willow is a principle factor driving differences in lignocellulosic sugar yield yet much of the physiology and development of this tissue is unknown. Traditional sectioning can be both difficult and destructive in woody tissue; however, technology such as three dimensional X-ray micro-computational tomography ( $\mu$ CT) scanning can be used to move biological researchers beyond traditional two dimensional assessment of tissue variation without having to destructively cut cells. This technology does not replace classical microscopic techniques but rather can be carefully integrated with traditional methods to improve exploration of the world of plant biology in three dimensions. The procedures below outline preparation of willow for 3D X-ray  $\mu$ CT and associated xylem staining and visualisation techniques, in particular secondary xylem programmed-cell-death (PCD) delay during gelatinous fibre (g-fibre) development. Many of the staining techniques here are transferable to other woody species such as poplar and Eucalyptus.

### Materials and Reagents

1. Glass slides
2. Cover slips
3. 48 well plates
4. Razor blades
5. Short rotation coppice (SRC) willow trees [cultivar Resolution-pedigree: {[*S. viminalis*. x (*S. viminalis*. x *S. schwerinii* SW930812)] x [*S. viminalis*. x (*S. viminalis*. x *S. schwerinii* 'Quest')]}]
6. 37% Formaldehyde (36.5-38% in H<sub>2</sub>O) (Sigma-Aldrich, catalog number: F8775)
7. Acetic acid ( $\geq 99.7\%$ ) (Sigma-Aldrich, catalog number: 320099)
8. Ethanol
9. FAA: 3.7% formaldehyde, 5% acetic acid and 47.5% ethanol
10. Oasis® floral foam (OASIS Floral Products) (<http://oasisfloralproducts.com/>)
11. Safranin-O solution (Sigma-Aldrich, catalog number: HT90432)
12. DPX Mountant for histology (Sigma-Aldrich, catalog number: 06522)

*Note: This product is named "Slide mounting medium (p-Xylene-bis(N-pyridinium bromide))" at Sigma-Aldrich.*


13. Distilled water
14.  $\beta$ -D-glucosyl Yariv reagent (Biosupplies Australia Pty Ltd., catalog number: 100-2)
15. HistoClear clearing agent (National Diagnostics, catalog number: HS-200)
16. Chlorazol black E/direct black 38 (Sigma-Aldrich, catalog number: C1144)

## **Equipment**

1. Reichert Sliding Microtome
2. Tweezers (for handling section)
3. Light microscope

## **Procedure**

X-ray  $\mu$ CT can also potentially be used for a broader range of crops, and a more diverse set of biological questions. In most cases, care should be taken to air-dry material prior to scanning; water saturation will reduce cell wall contrast while forced drying beyond fibre saturation point runs the risk of inducing cell collapse, altering the plant's tissue architecture. An example of this wider potential for the technology in plant research is provided in Figure 1 and Video 1.



**Figure 1. *Miscanthus* internode 3D render.** X-ray 3D Computational tomography example images of mature *Miscanthus x giganteus* stem.

**Video 1. 3D render of X-Ray  $\mu$ CT scans of willow (cultivar Resolution). 2 cm stem section (debarked). Vessel lumens are coloured in blue *in silico*.**




#### A. Stem sampling

The research where 3D X-ray  $\mu$ CT scanning was used to assess xylem fibre PCD (Brereton *et al.*, 2015) used short rotation coppice (SRC) willow trees (cultivar Resolution-pedigree: [*S. viminalis*. x (*S. viminalis*. x *S. schwerinii* SW930812)] x [*S. viminalis*. x (*S. viminalis*. x *S. schwerinii* 'Quest')]) grown in pots for 12 weeks. At week six the trees were tipped at a 45° angle to the horizontal and held in place for a further 6 weeks before harvesting (with the stem tied to a straight bamboo stem every two days to maintain the 45° angle). This alteration to the vector of gravitational stimulus elicits a unique change in tissue development known as the reaction wood response (Andersson-Gunneras *et al.*, 2003; Brereton *et al.*, 2011; Pitre *et al.*, 2007; Wardrop and Dadswell, 1955).

1. Harvest a sample whole stem segment of 6-8 cm from the midpoint of the stem.

*Note: This is deliberately chosen in preference to the traditional set breast height of 1.3 meters in order to allow developmentally comparable analysis between genotypes of different morphologies.*

2. Cut two small triangular notches into the wood 3-4 cm along the segment using a razor blade (making sure to reach through the bark to the secondary xylem), both notches should be at asymmetric positions from a transverse perspective (see Figure 2).



**Figure 2. Chlorazol black and Safranin O staining.** Classical sectioning and staining with light microscopy can be used to help navigate and orientate tissue architecture of the stem (Panel A). In this context, chlorazol black E (counter-stained with safranin O) is used to stain the g-layer of g-fibres in tension wood of a tipped stem sample of SRC willow (cultivar resolution)-this can then be aligned with images from the  $\mu$ CT scan (Panel B). Note, there is variation throughout the entirety of the 3D stem (longitudinal variation in tissue as well as transverse) so multiple notches in the stem biomass being scanned (one is visible here as a single triangle cut in panel A) are an important tool to help navigate the 3D volume data from the scan alongside natural asymmetry of the tissue. Scale bar = 4 mm.


3. Cut the segment in half to create two 3-4 cm segments. One 3-4 cm stem segment (stem A) will be used for X-ray  $\mu$ CT scanning while the other (stem B) will be used for associated staining and light or confocal microscopy. The notches will assist in navigation of the segments (allowing overlay and alignment) and can be removed *in silico* for 3D  $\mu$ CT render (stem A) or during sectioning (stem B).
4. To make sectioning, staining and X-ray  $\mu$ CT scan differentiation simpler, bark can be removed at this point from both stem segments by manual peeling or with a scalpel.
5. After harvesting stem segments should be “fixed” in FAA-formaldehyde, acetic acid, alcohol solution (3.7% formaldehyde, 5% acetic acid and 47.5% ethanol).

#### B. Drying and pre- $\mu$ CT scan prep

1. After fixation of Stem A for > 3 days, the segment should be removed from FAA solution and briefly washed in H<sub>2</sub>O 4 days prior to 3D X-ray  $\mu$ CT scanning.
2. Air-drying of stem A should be performed at room temperature (no greater than 30 °C) for over 4-5 days in a clean environment with relatively good air turnover. Forced drying at higher temperatures can pull H<sub>2</sub>O from the cell wall (beyond fibre saturation point) resulting in cell wall collapse. In very dry climates (low humidity) this can also occur during air-drying.



*Note: Air-drying of stem samples seems an essential step in creating contrast between the three elements within the biological system under scrutiny: the cell wall, the protoplast and void space. The important factor of air-dried samples being that “full” cells, tubular in shape, are easily distinguished from intracellular void space (which was not the case when “wet”, as intracellular void space is full of water (Figure 3).*



**Figure 3. Wet and air-dry biomass 2D computational tomography scan examples.**

3D X-Ray  $\mu$ CT models are rendered from >1,000 2D images, such as those from a tipped stem sample of SRC willow (cultivar resolution) shown above. A represents a 2D slice of a sample scanned “wet” whereas B represents a 2D slice of the same sample after having been air-dried for 3-4 days. Greater contrast between cell walls, cells with a protoplast (having not completed programmed cell death) and void volume is visible in the air-dried sample.

3. Once stem A is air-dried the sample can be secured in a low-density media, usually standard Oasis® floral foam (a phenolic based polymer, care should be taken if specifically studying cellulose, as such low densities are hard to differentiate) and positioned for scan with the structure/s of interest in the centre of the scan volume.

#### C. Infiltration

1. As X-ray absorption increases with atomic number of atoms within a voxel, stem samples can potentially be vacuum infiltrated with heavy metals before 3D X-ray  $\mu$ CT scanning. Such as with phosphotungstate (Staedler *et al.*, 2013).

#### D. Post- $\mu$ CT scan

1. Once stem A has been scanned (Brereton *et al.*, 2015) a destructive density measurement can be taken. Basic density (Biermann, 1996; Jourez *et al.*, 2001) is best suited for biomass density quantification as it incorporates tissue structural variation often meaningful in terms of physiology.

2. Vacuum infiltrate stem A with water before green volume is measured via water displacement (Biermann, 1996). Vacuum infiltration can be performed by pulling a vacuum over submerged wood samples for around 24 h depending on the material used. Wood oven dry weight is then measured after drying overnight at 105 °C.
3. Calculated as (Biermann, 1996; Jourez *et al.*, 2001): Basic density = oven dry mass/volume.

#### E. Sectioning and staining

Staining of the transverse sections from Stem B allows for different tissue types to be visualised and quantified within the stem. By clearly following each stem segment and navigating natural asymmetry (as well as notches cut into the stems), the *in silico* 3D X-ray  $\mu$ CT render, g-layers of g-fibres [chlorazol black E (Brereton *et al.*, 2015; Brereton *et al.*, 2011; Brereton *et al.*, 2012; Robards and Purvis, 1964)],  $\beta$ -galactosyl Yariv reagent (Tryfona *et al.*, 2012; Kitazawa *et al.*, 2013), secondary cell wall (safranin O) and remnants of cell contents (Coomassie and  $\beta$ -galactosyl Yariv) can be aligned and compared *in situ*. The *in situ*  $\beta$ -galactosyl Yariv reagent and Coomassie staining seem to be without published precedence but are simple and effective. Coomassie staining is particularly effective for macro navigation (across a whole transverse section) of cell content remnants using light microscopy, as classical viability staining can be ineffective in woody tissue (such as: Tunnel and NBT) due to cell disruption during the sectioning process. Staining duration and concentrations are only suggestions here as cell wall structure will vary depending on the plant studied and the growth conditions.


1. Once Stem B is secured in the Reichert sledge microtome, adjust angle of blade to the correct plane (this may require some trial and error depending on the material), dampen the sample with water and cut the thickness of < 30  $\mu$ m.
2. A 48 well plate is useful to perform multiple staining pipelines simultaneously on a common set of sections.
3. Using tweezers throughout, first place a section in H<sub>2</sub>O to hydrate for 2 min before using section for either chlorazol black E staining,  $\beta$ -galactosyl Yariv reagent *in situ* staining or Coomassie *in situ* staining (cell wall antibodies can also be integrated with analysis here).

#### **Chlorazol black E staining (example Figure 2)**

4. Counter-stain by placing sections in a well containing 1% safranin O (aq) for 2-3 min.
5. Wash in wells containing H<sub>2</sub>O x2 for 30 sec.
6. Place section in progressive dehydration series: 50% ethanol for 30 sec, 75% ethanol for 30 sec, 100% ethanol x2 for 1 min.
7. Place section in well containing 1% chlorazol black E (in methoxyethanol) for 4 min to stain g-layers.
8. Wash in well containing 100% ethanol x2 for 1 min.

### ***β*-galactosyl Yariv reagent staining (example Figure 4)**

4. Place section in 0.1 mg/ml dilute *β*-galactosyl Yariv reagent in 0.15 M NaCl and incubate at room temperature for 4 h on gently shaking platform.
5. Wash as necessary in wells containing H<sub>2</sub>O.
6. Move onto step E9.



**Figure 4. *In situ* *β*-galactosyl Yariv reagent staining.** *β*-galactosyl Yariv reagent has long been used to condense arabinogalactan proteins (AGPs). As the g-layer of g-fibres in tension wood are rich in Fasciclin-like arabinogalactan protein (Lafarguette *et al.*, 2004; Andersson-Gunneras *et al.*, 2006), *β*-galactosyl Yariv reagent works well as an *in situ* stain for g-layers in tension wood of SRC willow stem as well as a protoplast stain (live cells will

also include AGPs). Light microscopy of a tipped stem sample of SRC willow (cultivar resolution) shows A patterning of tensions wood g-layers and delayed PCD well aligned those of Figures 2. The same sample magnified using light microscopy to x63 illustrates: B. Tension wood (TW), C Opposite wood (OW) and D TW/OW interface show how clear the *in situ*  $\beta$ -galactosyl Yariv reagent staining can differentiate these cellular components. Scale bar = 4 mm.


#### ***Coomassie staining (example Figure 5)***

4. Place section in 1% Coomassie solution and incubate at room temperature for 4 h on gently shaking platform.
5. Wash as necessary in wells containing H<sub>2</sub>O.
6. Move onto step E9.

#### ***Slide construction***

9. Rinse staining section in Histoclear clearing agent for 1 min. At this stage apply DPX to slide to slightly set (1 min or so).
10. Place section on to the DPX and apply cover slip, lowering gently. Apply pressure to force out air bubbles, place some weight (a few grams) on top and leave to set for 24 h.





**Figure 5. *In situ* Coomassie staining.** A. Stitched 25 µm transverse sections stained with 1% Coomassie solution and counter stained with safranin O. B. Stitched 25 µm transverse sections without counterstain, remnants of cell contents is visible if present. While z-stacked images are difficult to assemble clearly using Coomassie staining, the process was invaluable to assess remnants of cell content using light microscopy (by eye, as one can move slightly through 3D space using light microscopes) and led to the identification of the same PCD patterning using 3D X-ray µCT. Scale bar = 4 mm.

### **Acknowledgments**

The author was financially supported by BioFuelNet Canada. X-ray µCT scanning was performed in collaboration with Farah Ahmed and Dan Sykes at the Natural History Museum with funding from BBSRC Sustainable Bioenergy Centre (BSBEC), working within the BSBEC BioMASS (<http://www.bsbec-biomass.org.uk/>) Programme. *In situ* staining techniques were developed in collaboration with Dr. Michael Jasmine Ray.

## References

1. Andersson-Gunneras, S., Hellgren, J. M., Bjorklund, S., Regan, S., Moritz, T. and Sundberg, B. (2003). [Asymmetric expression of a poplar ACC oxidase controls ethylene production during gravitational induction of tension wood.](#) *Plant J* 34(3): 339-349.
2. Andersson-Gunneras, S., Mellerowicz, E. J., Love, J., Segerman, B., Ohmiya, Y., Coutinho, P. M., Nilsson, P., Henrissat, B., Moritz, T. and Sundberg, B. (2006). [Biosynthesis of cellulose-enriched tension wood in Populus: global analysis of transcripts and metabolites identifies biochemical and developmental regulators in secondary wall biosynthesis.](#) *Plant J* 45(2): 144-165.
3. Biermann, C. J. (1996). Handbook of pulping and papermaking, 2nd edn. Academic Press.
4. Brereton, N. J., Ahmed, F., Sykes, D., Ray, M. J., Shield, I., Karp, A. and Murphy, R. J. (2015). [X-ray micro-computed tomography in willow reveals tissue patterning of reaction wood and delay in programmed cell death.](#) *BMC Plant Biol* 15: 83.
5. Brereton, N. J., Pitre, F. E., Ray, M. J., Karp, A. and Murphy, R. J. (2011). [Investigation of tension wood formation and 2,6-dichlorobenzonitrile application in short rotation coppice willow composition and enzymatic saccharification.](#) *Biotechnol Biofuels* 4: 13.
6. Brereton, N. J., Ray, M. J., Shield, I., Martin, P., Karp, A. and Murphy, R. J. (2012). [Reaction wood - a key cause of variation in cell wall recalcitrance in willow.](#) *Biotechnol Biofuels* 5(1): 83.
7. Jourez, B., Riboux, A. and Leclercq, A. (2001). [Comparison of basic density and longitudinal shrinkage in tension wood and opposite wood in young stems of Populus euramericana cv. Ghy when subjected to a gravitational stimulus.](#) *Canadian Journal of Forest Research-Revue Canadienne De Recherche Forestiere* 31, 1676-1683.
8. Kitazawa, K., Tryfona, T., Yoshimi, Y., Hayashi, Y., Kawauchi, S., Antonov, L., Tanaka, H., Takahashi, T., Kaneko, S., Dupree, P., Tsumuraya, Y. and Kotake, T. (2013). [beta-galactosyl Yariv reagent binds to the beta-1, 3-galactan of arabinogalactan proteins.](#) *Plant Physiol* 161(3): 1117-1126.
9. Lafarguette, F. et al. (2004). [Poplar genes encoding fasciclin-like arabinogalactan proteins are highly expressed in tension wood.](#) *New Phytologist* 164, 107-121.
10. Pitre, F. E., Cooke, J. E. K. and Mackay, J. J. (2007). [Short-term effects of nitrogen availability on wood formation and fibre properties in hybrid poplar.](#) *Trees-Structure and Function* 21, 249-259.
11. Robards, A. W. and Purvis, M. J. (1964). [Chlorazol black E as a stain for tension wood.](#) *Biotechnic and Histochemistry* 39, 309-315.



12. Staedler, Y. M., Masson, D. and Schonenberger, J. (2013). [Plant tissues in 3D via X-ray tomography: simple contrasting methods allow high resolution imaging.](#) *PLoS One* 8(9): e75295.
13. Timell, T. E. (1969). The chemical composition of tension wood. *Comptes Rendus Biologies* 72, 173-181.
14. Tryfona, T., Liang, H. C., Kotake, T., Tsumuraya, Y., Stephens, E. and Dupree, P. (2012). [Structural characterization of Arabidopsis leaf arabinogalactan polysaccharides.](#) *Plant Physiol* 160(2): 653-666.
15. Wardrop, A. and Dadswell, H. (1955). [The nature of reaction wood. IV. Variations in cell wall organization of tension wood fibres.](#) *Australian Journal of Botany* 3, 177-189.