

Target Gene Inactivation in Cyanobacterium *Anabaena* sp. PCC 7120

Kangming Chen¹, Huilan Zhu¹, Liping Gu¹, Shengni Tian^{1, 2} and Ruanbao Zhou^{1, *}

¹Department of Biology and Microbiology, South Dakota State University, Brookings, USA; ²College of Life Sciences, Anhui Agricultural University, Hefei, China

*For correspondence: Ruanbao.zhou@sdstate.edu

[Abstract] *Anabaena* sp. strain PCC 7120 has long served as a model organism for investigating N₂-fixation, photosynthesis, and various plant-type metabolic pathways and biofuel production, as well as cellular differentiation (Xu *et al.*, 2008, Halfmann *et al.*, 2014, Golden and Yoon, 2003). Since more than 30,000 sequenced bacterial genomes are currently available (Land *et al.*, 2015), specific gene inactivation and analyses of the corresponding mutant's phenotype have become powerful tools in elucidating the function of a target gene. Here we describe a protocol to inactivate a target gene in *Anabaena* sp. PCC 7120 using a single-crossover approach. This approach requires only one-step cloning of an internal fragment of a target gene into an integrative vector to produce a cargo plasmid. Upon a single crossover (homologous recombination) between the cargo plasmid and the *Anabaena* chromosome, the endogenous target gene is disrupted by generating 3'- and 5'-deleted fragments. This gene inactivating protocol is based on an integrative vector pZR606 (Chen *et al.*, 2015), which may be broadly applied to gene inactivation in other cyanobacterial species as well as other prokaryotic organisms.

Materials and Reagents

1. Petri dishes (Fisher Scientific, catalog number: FB0875713)
2. Glass flask (Fisher Scientific, catalog number: FB-500-50)
3. 1.5 ml centrifuge tubes (Fisher Scientific, catalog number: 02-682-550)
4. 15 ml conical centrifuge tubes (Fisher Scientific, catalog number: 05-527-90)
5. 50 ml conical centrifuge tubes (Fisher Scientific, catalog number: 12-565-270)
6. *Anabaena* sp. strain PCC 7120 (hereafter *Anabaena* 7120)
7. *Escherichia coli* (*E. coli*) NEB10 β [Δ (*ara-leu*) 7697 *araD*139 *fhuA* Δ *lacX*74 *galK*16 *galE*15 *e14- ϕ 80d**lacZ* Δ *M15* *recA*1 *relA*1 *endA*1 *nupG* *rpsL* (Str^R) *rph* *spoT*1 Δ (*mrr-hsdRMS-mcrBC*); New England BioLabs, catalog number: C3019H]
8. *E. coli* HB101 [F⁻*mcrB mrr hsdS*20(*r_B⁻*,*m_B⁻*) *recA*13, *leuB*6, *ara*-14, *proA*2, *lacY*1, *galK*2, *xyl*-5, *mtl*-1 *rpsL*20(Sm^R) *glnV*44 λ ⁻; Promega, catalog number: L2015]
9. Conjugal plasmid pRL443 and helper plasmid pRL623 (Elhai *et al.*, 1997)
10. pZR606, an integrative vector for *Anabaena* 7120 (Chen *et al.*, 2015). pZR606 is available upon request (GenBank, catalog number: KJ500179.1)
11. pZR670 (available upon request), a replicative vector in *Anabaena* 7120 (Xu *et al.*, 2005) [The

map of pZR670 is provided ([Figure S1](#)).]

12. LB broth (Sigma-Aldrich, catalog number: L3522-1KG)
13. LB agar, used for growing *E. coli* (MP, catalog number: 100262)
Note: This brand of agar does not work for cyanobacteria.
14. Agar, required for growing *Anabaena* 7120 and other cyanobacteria (Fisher Scientific, catalog number: A360-500)
15. Immobilon-NC transfer membrane (Millipore, catalog number: HATF08550)
16. Ampicillin sodium salt (Sigma-Aldrich, catalog number: A9518-25G)
17. Chloramphenicol (Fisher Scientific, catalog number: BP904-100)
18. Erythromycin (Sigma-Aldrich, catalog number: E6376-25G)
19. Kanamycin sulfate (Sigma-Aldrich, catalog number: K4000-25G)
20. Spectinomycin dihydrochloride pentahydrate (Sigma-Aldrich, catalog number: S9007-25G)
21. $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ (MP, catalog number: 194833)
22. $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ (Fisher Scientific, catalog number: BP510-500)
23. NaCl (Fisher Scientific, catalog number: S271-1)
24. K_2HPO_4 (Fisher Scientific, catalog number: BP363-500)
25. KNO_3 (Fisher Scientific, catalog number: BP368-500)
26. NaNO_3 (Fisher Scientific, catalog number: BP360-500)
27. $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$ (Fisher Scientific, catalog number: M87-100)
28. $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$ (99% purity) (Acros Organics, catalog number: 206371000)
29. $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ (Fisher Scientific, catalog number: Z76-500)
30. $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ (Fisher Scientific, catalog number: BP346-500)
31. H_3BO_3 (Fisher Scientific, catalog number: BP168-500)
32. NH_4VO_3 (Acros Organics, catalog number: 194910500)
33. $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$ (Fisher Scientific, catalog number: C371-100)
34. KOH (Fisher Scientific, catalog number: P250-500)
35. $\text{Na}_2\text{EDTA} \cdot 2\text{H}_2\text{O}$ (Fisher Scientific, catalog number: BP120-500)
36. $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ (Fisher Scientific, catalog number: I146-500)
37. Allen and Arnon medium plus nitrate: AA/8(N) (see Recipes)

Equipment

1. Centrifuges (Beckman, model: Allegra X-15R; Thermo Fisher Scientific, Sorvall™ Legend™ model: Micro17)
2. Innova-44R shaker equipped with continuous fluorescent light illumination (ca. 50-100 $\mu\text{E}/\text{m}^2 \text{ s}$) (Eppendorf, New Brunswick™ Innova®, model: 44R)
3. Cyanobacteria Culture Room (constant at 30 °C) equipped with a digital temperature controlling system and Sylvania fluorescent light bulbs (F40CWX, 40W/4100K, T12)

Note: Light shelves in cyanobacteria culture room are able to provide continuous light

illumination (ca. 50-150 $\mu\text{E}/\text{m}^2 \text{ s}$).

4. ESCO laminar flow cabinet (ESCO, Clean BenchAirstream® model: AHC-4B2)
5. Water bath sonicator (Bransonic® Ultrasonic Cleaner, model: 1510R-MT)
6. Microscope, Olympus upright compound microscope (Olympus, model: AX70; BX53)
7. PCR thermal cycler (Bio-Rad Laboratories, Touch™, model: C1000)
8. UV-Vis spectrophotometer (Thermo Fisher Scientific, model: GeneSyS 10S)

Procedure

A. Preparation of *E. coli* strains

1. Amplify the internal fragment (at least 500 bp) of the target gene from *Anabaena* chromosomal DNA.
2. Clone the internal fragment into the multiple cloning sites (*Bgl*II-*Not*I-*Age*I-*Spe*I-*Apa*I-*Sma*I/*Xma*I) of pZR606 (Km^R/Sp^R) to produce the cargo plasmid [see details in (Chen *et al.*, 2015)]. The map for multiple cloning sites in pZR606 is provided ([Figure S2](#)). The cargo plasmid must be transformed to *E. coli* NEB10 β or *E. coli* DH10B.
3. Grow *E. coli* HB101 harboring pRL443 (Ap^R) and pRL623 (Cm^R) in 2 ml LB broth with 100 $\mu\text{g}/\text{ml}$ ampicillin and 25 $\mu\text{g}/\text{ml}$ chloramphenicol overnight (~16 h) at 37 °C with 200 rpm shaking. Meanwhile, grow NEB10 β bearing the cargo plasmid (Km^R/Sp^R) in 2 ml LB broth with 50 $\mu\text{g}/\text{ml}$ kanamycin overnight at 37 °C with 200 rpm shaking.
4. Transfer 100 μl of the above overnight *E. coli* cultures into 5 ml fresh LB with appropriate antibiotics respectively; continue to grow for ~3 h ($\text{OD}_{600} \sim 0.5$).
5. Harvest the 5 ml culture by centrifugation at 4,000 $\times g$ in Allegra X-15R for 10 min at 25 °C (room temperature).
6. Wash the cell pellets three times with 1 ml LB to remove antibiotics completely, then add 200 μl LB to re-suspend the pellets, respectively.
7. For mating experiment, mix 100 μl NEB10 β bearing cargo plasmid with 100 μl HB101 harboring pRL443 and pRL623 and incubate at room temperature for 30 min (experimental group). For a negative control, mix 100 μl LB with 100 μl HB101 harboring pRL443 and pRL623 and incubate at room temperature for 30 min (control group).

B. Preparation of *Anabaena* 7120 culture

8. Grow *Anabaena* 7120 in 30 ml AA/8(N) medium (Allen & Arnon, 1955a, Allen & Arnon, 1955b) or BG11 medium (Rippka *et al.*, 1979) for ca. 5 days until it reaches early exponential stage ($\text{OD}_{700} \sim 0.5$). The *Anabaena* 7120 is grown under continuous light illumination (ca. 50 $\mu\text{E}/\text{m}^2 \text{ s}$) at 30 °C and shaken at 120 rpm in a temperature controlled Innova-44R lighted incubator (New Brunswick Scientific) or in the Cyanobacteria Culture Room with continuous light illumination (ca. 50 $\mu\text{E}/\text{m}^2 \text{ s}$) at 30 °C and shaken at 120 rpm.
9. Harvest the culture by centrifugation at 4,000 $\times g$ with Allegra X-15R for 10 min at 25 °C.

10. Re-suspend the cell pellet with 2 ml AA/8(N) in a 25 ml glass flask, break filaments into an average 3-5 cell lengths (confirmed microscopically) by sonicating the cultures for 60-120 sec using water bath sonicator ([Figure S3](#)). Set up the sonicator under standard model and add water to the operating level line during the sonication process.
 11. Transfer the sonicated culture into a 15 ml centrifuge tube and centrifuge in Allegra X-15R at 4,000 x g for 10 min at 25 °C. Re-suspend the cell pellet with 1 ml AA/8(N), transfer the cells into a 1.5 ml Eppendorf tube, and then centrifuge at 6,000 x g for 1 min (Thermo Fisher Scientific, Micro17) at 25 °C. Re-suspend the pellet with 400 µl AA/8(N); then divide cells equally into two tubes.
- C. Conjugal transformation of a cargo plasmid into *Anabaena* 7120
12. Mix 200 µl of *Anabaena* 7120 resuspension (step 11) with 200 µl of the mated *E. coli* mixture (experimental group, step 7), or the control mixture (control group, step 7), respectively, and incubate at room temperature for another 30 min.
 13. Spread the conjugal mixtures onto the autoclaved Immobilon-NC transfer membrane atop AA (N) agar containing 5% LB (v/v) without antibiotic and incubate at 30 °C, with continuous light illumination (ca. 50 µE/m² s) for 24 h.
 14. Transfer the membrane onto AA (N) agar plate with 10 µg/ml spectinomycin (Sp10), and incubate under the same growth conditions until spectinomycin-resistant colonies gradually appear during 10-15 days as shown in Figure 1 (top and middle panels).




Figure 1. Screen for the potential knockout mutants. The conjugal transfer mixture with the cargo plasmid was spread onto Immobilon-NC membrane atop AA (N) agar containing Sp10 for 0 (*top panel*) and 19 days (*middle panel*). The bottom panel shows the above individual colonies re-streaked onto a new AA (N) agar plate containing Sp10 for segregation. As a negative control, the conjugal transfer mixture without cargo plasmid (control) showed no individual colonies on AA (N) agar plate containing Sp10 (data not shown).

D. Verification of the knockout mutants

15. Re-streak at least ten potential knockout mutant colonies onto a fresh AA (N) agar plate containing Sp10, allowing complete segregation of mutant chromosomes (Figure 1, bottom panel) since *Anabaena* 7120 has multiple copies of chromosomes per cell (Hu *et al.*, 2007).
16. Screen for the positive knockout mutant colonies by colony PCR. Three primers were designed to verify the target gene knockout mutants. Targeted gene open reading frame (ORF) primer pair [forward primer (FP) and reverse primer (RP)] was used to verify the wild-type gene; while forward primer FP paired with a vector-specific primer ZR90 (AAGTTCTTCTCCTTTGCTAGC) (Chen *et al.*, 2015) was used to distinguish the positive knockout mutants from WT. A standard PCR program was used with the exception that the chromosomal DNA template was prepared

- by preheating cyanobacterial cell suspension in 6 μ l ddH₂O at 95 °C for 10 min.
17. Completely suspend PCR positive mutant colonies in 100 μ l AA/8(N) medium to make serial dilutions (1x, 10x and 100x). Then re-streak each dilution onto AA (N) agar plate containing Sp10 for single colonies to separate *E. coli* contamination from *Anabaena*. These single *Anabaena* knockout mutant colonies were confirmed again with the colony PCR described above.
 18. Repeat steps 15 to 17 until no wild-type target gene was detected in the knockout mutant colony PCR. An example of colony PCR verification for the inactivated mutants of *all4160* in *Anabaena* 7120 is shown in Figure 2.




Figure 2. Colony PCR verification for the single crossover inactivated *all4160* mutant in *Anabaena* 7120. The target gene specific F₄₁₆₀ + R₄₁₆₀ and F₄₁₆₀ + ZR90 primer pairs were used for PCR amplification from the knockout mutant colonies of *all4160* (lanes 2-3 and 5-6, respectively) or from the WT 7120 (lanes 4 and 7) as controls. The PCR product sizes amplified with F₄₁₆₀ + R₄₁₆₀ and F₄₁₆₀ + ZR90 primer pairs are expected to be 1,056 bp (lane 7) and 910 bp (lanes 2-3), respectively. The sizes of DNA ladder (lane 1) are also indicated. Lanes 2 & 3 or lanes 5 & 6 were two independent knockout mutants. Lanes 5 & 6 show that the target gene was successfully knocked out, the mutants were completely segregated. F₄₁₆₀: ggatccatgTACATGGCAACCAAAGTG, R₄₁₆₀: gtcgacATAAGCGCCACTATTTCTATTAAA (uppercase letters represent nucleotides corresponding to *all4160*).

E. Complementation experiment

19. A complementation experiment is required to confirm that the target gene, rather than a downstream gene, is responsible for the mutant phenotype, as in some cases, a knockout mutant's phenotype might be due to a polar effect on the expression of a downstream gene. Complementation by a single gene was tested by sub-cloning the entirety of the target gene-coding region into a shuttle vector pZR670 (Xu *et al.*, 2015) or pRL2833a (Wolk *et al.*, 2007), which contains both the chloramphenicol (Cm^R) and the erythromycin (Em^R) resistance genes. To construct the complementing plasmid, the coding sequence of the target gene is cloned into the multiple cloning sites (MCS) (*Nsi*I/*Nde*I/*Avr*II/*Aat*II/*Xho*I/*Asc*I/*Nae*I/*Bam*HI) of pZR670 (Figure 3), in-frame with the glnA start codon. For example, coding sequencing of *alr4853* was in-frame cloned into *Nde*I-*Xho*I digested pZR670, producing pZR1617 (Xu *et al.*, 2015).




Figure 3. The multiple cloning sites of a complementing vector pZR670. P_{glnA}: promoter region of *glnA* (*alr2328*), RBS: Ribosomal binding site, 5'-glnA: N-terminal portion (MTTPQEVLKRI) of *Anabaena* GlnA (glutamine synthetase), MCS (nucleotides highlighted in red): multiple cloning sites. The target gene-coding sequence (ORF) can be in-frame inserted between *NdeI* and one of *AvrII/AatII/XhoI/AscI/NaeI/BamHI* digested pZR670 (Xu et al., 2015).

20. A similar approach described in Part C is used to transfer the complementing plasmid into the knockout mutant. *E. coli* and *Anabaena* mutant culture were prepared as described before (steps 3 to 11), except that the cargo plasmid was substituted by the complementing plasmid. Then dilute the mutant *Anabaena* culture by 5x, 10x and 100x with AA/8(N) in 1.5 ml centrifuge tubes. Spot 5 µl of the diluted cultures onto the Immobilon-NC membrane atop AA(N) agar containing 5% LB without antibiotic. Make 3 replicate spots for each sample. Air-dry the membrane for a few minutes; then spot 5 µl of the mated *E. coli* mixture atop the *Anabaena* spots. Air-dry for a few minutes, and incubate the plate under continuous light illumination (ca. 50 µE/m² s) at 30 °C for 24 h.
21. Transfer the membrane onto AA (N) agar plate containing both Sp10 and 10 µg/ml of erythromycin, and incubate under the same growth conditions. The Sp-Em-resistant colonies will grow in 10-15 days, similarly as seen in Figure 1.
22. Then verify the Sp-Em-resistant colonies that contain the complementing plasmid by colony PCR with specific primers targeting the Em resistance gene.
23. Examine the phenotype for at least three verified colonies to see if the complementing plasmid is capable of complementing the single crossover generated mutant. If the complemented mutant displays a restored wild-type phenotype, in which the target gene has been successfully inactivated by a single crossover approach, you can conclude that the target gene is responsible for the mutant's phenotype, which hints the function of the target gene.
24. If the complementing plasmid fails to restore the mutant's phenotype to its wild-type's phenotype in step 22, this indicates that the mutant's phenotype might be caused by a polar effect on a downstream gene. No functional clue could be made for the target gene under this condition although it has been successfully inactivated.

Notes

1. *Anabaena* sp. strain PCC 7120 can also be grown in BG11 medium (<http://www-cyanosite.bio.purdue.edu/media/table/BG11.html>) (Rippka et al., 1979).
2. After transferring the membrane onto AA (N) agar plate with appropriate antibiotics, keep tracking the membrane color change. The wild-type *Anabaena* 7120 should be killed and the color of the membrane should gradually turn to yellow.
3. Change the AA (N) agar plate with appropriate antibiotics every week to maintain effective antibiotics. If the membrane in the control plate (*Anabaena* 7120 + HB101 harboring pRL443 and pRL623) does not turn to yellow in one week, transfer all the membranes onto new AA(N) agar plate containing appropriate antibiotics immediately.

Recipes

1. Allen and Arnon medium plus nitrate: AA/8(N)
Anabaena 7120 is grown in a modified Allen-Arnon medium (Allen & Arnon, 1955a). AA liquid medium was diluted 8 times with sterilized distilled H₂O to make AA/8. For making AA/8(N), a final concentration of 2.5 mM nitrate was added to AA/8 medium from KNO₃/NaNO₃ stock solution. For preparation of AA(N) agar plates, a final concentration of 2.5 mM nitrate and 1% agar were added to AA medium.
 - a. AA medium preparation
Solution A 25.0 ml
Solution B 6.25 ml
Agar (if needed) 10.0 g
Distilled water 969.0 ml
Aliquot 200 ml into five bottles, and then autoclave them at 121 °C for 20 min
 - b. Solution A (A & A minus phosphate stock solution)
4% MgSO₄·7H₂O 500.0 ml (final concentration 1%, m/m)
1.2% CaCl₂·2H₂O 500.0 ml (final concentration 0.3%, m/m)
3.8% NaCl 500.0 ml (final concentration 0.95%, m/m)
Microelements stock solution 500.0 ml (final concentration of microelements in solution A is four times diluted microelements stock solution)
 - c. Solution B (K₂HPO₄ stock solution)
K₂HPO₄ 21.40 g
Distilled water 500.0 ml
 - d. KNO₃/NaNO₃ stock solution (500 mM)
KNO₃ 25.276 g
NaNO₃ 21.249 g
Add distilled H₂O to a final volume of 500 ml

Stock solution is autoclaved and stored at 4 °C

e. Microelements stock solution

A & A FeEDTA solution 160.0 ml

MnCl₂·4H₂O 360 mg

Na₂MoO₄·2H₂O 61.1 mg

ZnSO₄·7H₂O 44.0 mg

CuSO₄·5H₂O 15.8 mg

H₃BO₃ 572.0 mg

NH₄VO₃ 4.6 mg

CoCl₂·6H₂O 8.0 mg

Distilled water 1,090.0 ml

f. A & A FeEDTA stock solution

Dissolve 5.2 g KOH in 186 ml distilled water, add 20.4 g Na₂EDTA·2H₂O

Dissolve 13.7 g FeSO₄·7H₂O in 364 ml distilled water

Mix the above two solutions, then bubbling millipore-filtered air through solution until color changes. The final pH of FeEDTA solution is approximately 7.5.

Note: All stock solutions are recommended to store at 4 °C.

Acknowledgments

The protocol is based on the publications “Conjugal transfer of DNA to cyanobacteria” (Elhai and Wolk, 1988); “Simultaneous gene inactivation and promoter reporting in cyanobacteria” (Chen *et al.*, 2015), and “Characterization of five putative aspartate aminotransferase genes in the N₂-fixing heterocystous cyanobacterium *Anabaena* sp. strain PCC 7120” (Xu *et al.*, 2015). The authors would like to thank Jaimie Gibbons for her critical reading of the manuscript. This work was partially supported by the NSF, Energy for Sustainability Grant CBET1133951 (to R. Z.), and by the USDA-NIFA grant 11665597 (to R. Z.).

References

1. Allen, M. B. and Arnon, D. I. (1955). [Studies on nitrogen-fixing blue-green algae. I. growth and nitrogen fixation by *Anabaena cylindrica* Lemm.](#) *Plant Physiol* 30(4): 366-372.
2. Chen, K., Xu, X., Gu, L., Hildreth, M. and Zhou, R. (2015). [Simultaneous gene inactivation and promoter reporting in cyanobacteria.](#) *Appl Microbiol Biotechnol* 99(4): 1779-1793.
3. Elhai, J., Vepritskiy, A., Muro-Pastor, A. M., Flores, E. and Wolk, C. P. (1997). [Reduction of conjugal transfer efficiency by three restriction activities of *Anabaena* sp. strain PCC 7120.](#) *J Bacteriol* 179(6): 1998-2005.
4. Elhai, J. and Wolk, C. P. (1988). [Conjugal transfer of DNA to cyanobacteria.](#) *Methods Enzymol* 167: 747-754.

5. Golden, J. W. and Yoon, H. S. (2003). [Heterocyst development in *Anabaena*](#). *Curr Opin Microbiol* 6(6): 557-563.
6. Halfmann, C., Gu, L. and Zhou, R. (2014). [Engineering cyanobacteria for the production of a cyclic hydrocarbon fuel from CO₂ and H₂O](#). *Green Chemistry* 16: 3175-3185.
7. Hu, B., Yang, G., Zhao, W., Zhang, Y. and Zhao, J. (2007). [MreB is important for cell shape but not for chromosome segregation of the filamentous cyanobacterium *Anabaena* sp. PCC 7120](#). *Mol Microbiol* 63(6): 1640-1652.
8. Land, M., Hauser, L., Jun, S. R., Nookaew, I., Leuze, M. R., Ahn, T. H., Karpinets, T., Lund, O., Kora, G., Wassenaar, T., Poudel, S. and Ussery, D. W. (2015). [Insights from 20 years of bacterial genome sequencing](#). *Funct Integr Genomics* 15(2): 141-161.
9. Rippka, R., Deruelles, J., Waterbury, J. B., Herdman, M. and Stanier, R. Y. (1979). [Generic assignments, strain histories and properties of pure cultures of cyanobacteria](#). *Microbiology* 111: 1-61.
10. Wolk, C. P., Fan, Q., Zhou, R., Huang, G., Lechno-Yossef, S., Kuritz, T. and Wojciuch, E. (2007). [Paired cloning vectors for complementation of mutations in the cyanobacterium *Anabaena* sp. strain PCC 7120](#). *Arch Microbiol* 188(6): 551-563.
11. Xu, X., Elhai, J. and Wolk, C. (2008). [Transcriptional and developmental responses by *Anabaena* to deprivation of fixed nitrogen](#). Caister Academic Press, Norfolk, United Kingdom.
12. Xu, X., Gu, L., He, P. and Zhou, R. (2015). [Characterization of five putative aspartate aminotransferase genes in the N₂-fixing heterocystous cyanobacterium *Anabaena* sp. strain PCC 7120](#). *Microbiology* 161(6): 1219-1230.