

www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

Extraction and Assays of ADP-glucose Pyrophosphorylase, Soluble Starch Synthase and Granule Bound Starch Synthase from Wheat (*Triticum aestivum* L.) Grains

Konstantin Kulichikhin^{1, 2}, Shalini Mukherjee¹ and Belay T. Ayele^{1, *}

¹Department of Plant Science, University of Manitoba, Winnipeg, Manitoba, Canada; ²Laboratory of Genomic and Proteomic Researches, The Institute of Translational Biomedicine, Saint Petersburg University, Saint-Petersburg, Russia

*For correspondence: <u>belay.ayele@umanitoba.ca</u>

[Abstract] Starch biosynthesis in plants involves a network of enzymes of which adenosine-5'diphosphoglucose (ADP-glucose) pyrophosphorylase (AGPase, E.C. 2.7.7.27), and soluble and granule bound starch synthases (SSS and GBSS, E.C. 2.4.1.21) play central roles. Here, we outline the protocol for extraction and assay of these enzymes in developing grains of wheat (*Triticum aestivum* L.). The principle of the assays outlined is based on a coupling enzymatic reactions where the product of the initial reaction is used as a substrate for subsequent reactions in order to generate NADPH, which can be measured easily by spectrophotometer. This protocol does not need expensive labelled chemicals and can be carried out using equipment commonly found in a biochemical laboratory. We applied this protocol to study the dynamics of AGPase, SSS and GBSS activity in developing wheat grains at different time points after anthesis.

[Background] Starch is a carbohydrate polymer made up of amylose, a linear glucan polymer composed of α -1,4-linked glucose molecules, and amylopectin, another glucan polymer composed of α -1,4-linked glucose molecules branched by α -1,6-glycosidic bonds. The enzyme adenosine-5'-diphosphoglucose (ADP-glucose) pyrophosphorylase (AGPase, E.C. 2.7.7.27) catalyzes the first committed step of starch synthesis, converting glucose-1-phosphate and ATP to ADP-glucose and inorganic pyrophosphate (PPi). ADP-glucose is subsequently used by soluble starch synthases (SSS) and granule bound starch synthases (GBSS) (E.C. 2.4.1.21), and starch branching enzymes to elongate and branch the glucan chains of the starch granule.

Initially, AGPase and starch synthase assays were carried out using ¹⁴C- and ³²P-labelled ADPglucose (Ghosh and Price, 1966; Vos-Scheperkeuter *et al.*, 1986), which requires the use of expensive chemicals as well as specialized equipment to work with labelled compounds. Here, we outline the method adopted and applied for extracting and assaying AGPase, SSS and GBSS activity in developing wheat (*Triticum aestivum* L.) grains (Mukherjee *et al.*, 2015). Our protocol is based on the methods reported previously by Nakamura *et al.* (1989) for measuring AGPase and SSS activities in the endosperm of developing rice grains, and by Schaffer and Petreikov (1997) for measuring the activities of SSS and GBSS in tomato fruits. The method is based on coupling enzymatic reactions where the product of the initial reaction is used as a substrate for subsequent reactions in order to generate NADPH, which can be easily measured by spectrophotometers. The protocol can be carried out using the


www.bio-protocol.org/e1929


Vol 6, Iss 18, Sep 20, 2016

equipment commonly found in biochemical laboratories and does not require the use of labeled compounds.

Because this protocol is an adoption of the methods developed for different tissues of other plant species (endosperm of rice grain and tomato fruit), the reaction mixture composition for all enzymes was optimized for reaction buffer pH and substrate concentration so that enzyme activity was within the linear phase with respect to incubation time and protein concentration. The amounts of enzyme preparation added (PGM, pyruvate kinase, hexokinase, G6PDH) have been adjusted to achieve completion of coupling reaction in expected time frame. This protocol can be used for studying the activities of AGPase, SSS and GBSS in other tissues of wheat as well as in different tissues of other plant species; however, optimization of reaction buffer pH and substrate concentrations is required.

The principle of AGPase assay is presented in Figure 1. AGPase present in plant tissue extract catalyzes the conversion of ADP-glucose and PPi into glucose-1-phosphate and ATP. After the first reaction is stopped by boiling, phosphoglucomutase (PGM) is added to the reaction mixture; PGM converts quantitatively glucose-1-phosphate into glucose-6-phosphate. Subsequently, glucose-6-phosphate dehydrogenase (G6PDH) is added; G6PDH converts glucose-6-phosphate in the presence of NADP into 6-phosphogluconic acid. At the same time, NADP is converted into NADPH and the quantity of NADPH formed is equivalent to the quantity of glucose-6-phosphate oxidized. By measuring NADPH concentration spectrophotometrically at 340 nm, we can quantify the amount of ADP-glucose degraded by AGPase activity.





The principle of SSS/GBSS assay is presented in Figure 2. Starch synthases present in plant tissue extract (SSS) or in starch granule suspension (GBSS) catalyze the conversion of ADP-glucose into ADP coupled with the elongation of amylopectin primer for one glucose residue. Pyruvate kinase added to the reaction mixture after the first reaction is stopped by boiling converts quantitatively ADP into ATP in the presence of PEP. Subsequently, addition of glucose together with hexokinase leads to the conversion of ATP into glucose-6-phosphate. In the next step of the assay, G6PDH converts glucose-6-phosphate in the presence of NADP into 6-phosphogluconic acid; NADP is converted into NADPH at the same time and the quantity of NADPH formed is equivalent to the quantity of glucose-6-phosphate



oxidized. By measuring NADPH concentration spectrophotometrically at 340 nm, we can quantify the amount of ADP-glucose (at the first step of the analysis) degraded by SSS/GBSS activity.




Figure 2. Principle for soluble and granule bound starch synthases (SSS and GBSS) assay

Materials and Reagents

- A. Chemicals
 - 1. Common chemicals, buffers and kits
 - a. Protein assay kit (Bio-Rad Laboratories, catalog number: 5000002)
 - b. 4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid (HEPES) (Thermo Fisher Scientific, Fisher Scientific, catalog number: BP310)
 - c. Sodium hydroxide (NaOH) (Thermo Fisher Scientific, Fisher Scientific, catalog number: BP359)
 - d. Magnesium chloride hexahydrate (MgCl₂·6H₂O) (Sigma-Aldrich, catalog number: M2670)
 - e. Ethylenediaminetetraacetic acid (EDTA), disodium salt dihydrate (Na₂EDTA·2H₂O) (Sigma-Aldrich, catalog number: E1644)
 - f. DL-dithiothreitol (DTT) (Sigma-Aldrich, catalog number: D0632)
 - g. Polyvinylpyrrolidone (average molecular weight 10,000) (Sigma-Aldrich, catalog number: PVP10)
 - h. Potassium chloride (KCI) (Thermo Fisher Scientific, Fisher Scientific, catalog number: BP366)
 - i. Tetrasodium pyrophosphate (PPi) decahydrate (Na₄P₂O₇·10H₂O) (Sigma-Aldrich, catalog number: S6422)
 - 2. Substrates and cofactors
 - a. Adenosine-5'-diphosphoglucose (ADP-glucose), disodium salt (Sigma-Aldrich, catalog number: A0627)
 - b. Amylopectin from maize (Sigma-Aldrich, catalog number: 10120)

Copyright © 2016 The Authors; exclusive licensee Bio-protocol LLC.



www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

- c. β-Nicotinamide adenine dinucleotide phosphate (NADP) sodium salt hydrate (Sigma-Aldrich, catalog number: N0505)
- d. Phospho(enol)pyruvic acid (PEP) monosodium salt hydrate (Sigma-Aldrich, catalog number: P0564)
- 3. Enzymes
 - a. Hexokinase (E.C. 2.7.1.1) from *Saccharomyces cerevisiae* (Sigma-Aldrich, catalog number: H5000)
 - b. Glucose-6-phosphate dehydrogenase (G6PDH, E.C. 1.1.1.49) from *Saccharomyces cerevisiae* (Sigma-Aldrich, catalog number: G7877)
 - c. Phosphoglucomutase (PGM, E.C. 5.4.2.2) from rabbit muscle (Sigma-Aldrich, catalog number: P3397)
 - d. Pyruvate kinase (E. C. 2.7.1.40) from rabbit muscle (Sigma-Aldrich, catalog number: P1506)
- 4. Stock solutions (see Recipes)
 - a. 20 mM ADP-glucose
 - b. 10 mM NADP
 - c. 20 mM PEP
- B. Enzymatic activities reaction mixtures
 - 1. Extraction buffer (see Recipes, section 1)
 - 2. AGPase activity reaction mixture (see Recipes, section 3a)
 - 3. Starch synthase activity reaction mixture (see Recipes, section 3b)
 - 4. Solution 1 for starch synthase activity assay (Pyruvate kinase reaction mixture) (see Recipes, section 3c)
 - 5. Solution 2 for starch synthase activity assay (Glucose-6-phosphate dehydrogenase reaction mixture) (see Recipes, section 3d)
- C. Plant material
 - 1. Grains of common wheat (Triticum aestivum L.)
- D. Other materials
 - 1. Eppendorf tubes (1.5 and 2.0 ml)
 - 2. Falcon tubes (15 ml)
 - 3. Liquid nitrogen
 - Microcentrifuge tube locker (for example, Sigma-Aldrich, catalog number: Z708372, or Sorenson Bioscience MCT LidLock[™], catalog number: 11870)
 Note: Product Z708372 has been discontinued.
 - 5. Miracloth (Merck Millipore, catalog number: 475855)



Vol 6, Iss 18, Sep 20, 2016

Equipment

- 1. Porcelain mortar and pestle
- 2. Analytical balance (capacity 100 g or higher, resolution 0.001 g or better)
- 3. Digital heat block or water bath suitable for 30 °C and digital heat block or water bath suitable for 100 °C

Note: For some steps of the analysis, immediate transfer from 30 °C to 100 °C is necessary. This is why two heating devices are necessary simultaneously for the protocol.

- 4. Refrigerating ultracentrifuge with appropriate rotor suitable for spinning of 15 ml Falcon tube and 1.5-2 ml Eppendorf tubes at up to $12,000 \times g$
- Ultraviolet (UV)/Visible spectrophotometer equipped with heated cell holder suitable for maintaining cuvette temperature at 30 °C. We used GE Healthcare/Amersham Bioscience Ultrospec 3100 Pro UV/Visible spectrophotometer (Biochrom Ltd, Cambridge, UK).



www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

Procedure

Summary of the enzyme extraction procedure is shown in a flowchart (Figure 3).




Figure 3. Flowchart for extraction of AGPase, SSS and GBSS from developing wheat seeds. Please, refer to the main text of the protocol for specific conditions and buffers/media composition required for each step.

www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

A. Wheat grain harvesting and storage before analysis

Developing grains of wheat were harvested from primary and secondary tillers at different time points after anthesis (Figure 4), in liquid nitrogen and stored at -80 °C until further use. Care has to be taken not to thaw the frozen grains till enzyme assay as thawing the grains would seriously affect the result of the analysis.



Figure 4. Developing seeds of wheat cv. Yecora Rojo at 4, 8, 16 and 25 days after anthesis

- B. Extraction
 - 1. Wheat grains were ground to a fine powder using liquid nitrogen in a pre-cooled mortar and pestle and ground tissue was transferred into empty pre-weighed Falcon tube pre-cooled with liquid nitrogen. The weight of the ground tissue was determined. This was performed by first weighing an empty Falcon tube, and then weighing it again filled with the ground tissue. The difference of the two readings equals to the net weight of the ground tissue. This value was recorded and used for further calculations (*m*, see Calculations, Equations 1, 2, 4, 5, 7).
 - 2. The extraction buffer was prepared and kept on ice (see Recipes, section 1).
 - 3. 1 ml of ice-cold extraction buffer was added to each 100 mg of the ground tissue and then homogenized.
 - 4. The homogenate was filtered through two layers of Miracloth moistened with extraction buffer and centrifuged at 12,000 *x g* for 10 min at 4 °C, and the supernatant containing the crude enzyme was transferred to a fresh tube and stored on ice.
 - 5. The yellow layer of storage protein that forms over the pellet of starch granules was scraped off and discarded (Figure 5).

bio-protocol www.bio-protocol.org/e1929 Vol 6, Iss 18, Sep 20, 2016 -AGPase and SSS assay Centrifugation 10000 g, 10 min Enzyme Supernatant extract ne extract) Storage laver Removal of storage Strarch protein layer by spatula Homogenate granules GBSS laver assay Starch granules Pellet (storage proteins layer above starch granules layer)

Figure 5. Fractioning of plant tissue homogenate. After centrifugation, the pellet consists of two layers; the upper layer formed by storage proteins and the lower layer by starch granules. The supernatant is collected into a separate tube and is used for AGPase and SSS assay. The upper layer of the pellet (storage proteins) is removed by spatula and discarded. The lower layer (starch granules) is used for GBSS assay after washing.

6. The pellet of starch granules was washed with 2-4 ml* of extraction buffer twice by centrifugation to remove all traces of SSS. Supernatant obtained after first washing was combined with initial one (step B4) and the pooled supernatant was used for the AGPase and SSS activity assays. The volume of the pooled supernatant was recorded for further calculations (*V*₁, see Calculations, Equations 1, 2, 4, 5). Supernatant after second washing was discarded**.

*Note: The volume of extraction buffer used for washing should be sufficient to wash the starch granule pellet, but not in excess to avoid dilution of the extract. Use of extraction buffer volume that is equal to 20-25% of the volume used in step B3 is a good practice.

**Note: The aim of starch granule pellet washing is to eliminate SSS contamination of the GBSS fraction. The amount of SSS in the first wash may significantly affect the final result of the SSS activity assay. This is why the supernatant from the first wash is combined with the main extract. The second wash may contain some SSS activity, but this amount does not significantly affect the final result of the assay and it can be discarded to avoid extra dilution of the main extract.

The pellet of starch granules was re-suspended in 0.5-2.0 ml* of extraction buffer and stored on ice until used for the GBSS activity assay. The volume of suspension was recorded and used for further calculations of GBSS activity (*V*₆, see Calculations, Equation 7).
 *Note: The volume of buffer used to re-suspend starch granules depends on the amount of

granules collected. For larger quantities of granules, use a larger volume of re-suspending buffer, but do not exceed the volume used in step B6.

Copyright © 2016 The Authors; exclusive licensee Bio-protocol LLC.



www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

C. Enzyme activity assay

1. AGPase activity assay

Summary of the procedures for AGPase activity assay is shown in a flowchart (Figure 6).




Figure 6. Flowchart of AGPase assay. Please, refer to the main text of the protocol for specific conditions and buffers/media composition needed for each step of the assay.

Copyright © 2016 The Authors; exclusive licensee Bio-protocol LLC.

www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

- a. Two tubes with AGPase reaction mixture, each containing 50 mM HEPES-NaOH (pH 7.5), 1.2 mM ADP-glucose*, 5 mM PPi, 6 mM MgCl₂, 3 mM DTT and 50 µl of enzyme extract (V₂, see Calculations, Equations 1, 3) in a total volume of 0.45 ml were prepared.
 *Note: The reaction mixture is prepared first by mixing all components other than ADP-glucose in a volume of 423 µl as shown in Table 1 (see Recipes, section 3a).
- b. The tubes were pre-incubated at 30 °C for 2 min prior to starting the reaction by adding 27 μ I of 20 mM ADP-glucose stock solution.
- c. The reaction in the first tube (reference tube) was stopped immediately after the addition of ADP-glucose by incubating the tube in water bath at 100 °C for 1 min, whereas the second tube (experimental tube) was incubated at 30 °C for 15 min (*t*, reaction time, see Calculations, Equations 1, 3) followed by incubation of the tube at 100 °C for 1 min*. Subsequently (steps C1d-f), the reactions in the experimental and reference tubes were processed simultaneously in a similar manner.

*Note: Use special microcentrifuge tube lockers (see Materials and Reagents, section D4) to prevent tube opening during boiling. Tube opening at this stage may lead to partial evaporation of water and uncontrolled decrease of reaction mixture volume that may seriously affect the accuracy of the measurement.

- d. The reaction mixtures were allowed to cool down to room temperature followed by centrifuging at $10,000 \times g$ at room temperature for 10 min.
- e. The supernatant (0.45 ml) was transferred to a new tube, and 0.1 ml of 6 mM NADP and 0.3 ml of 50 mM HEPES-NaOH (pH 7.5) were added to it. The total volume of the final reaction mix (0.85 ml, *V*₃, see Calculations, Equations 1, 3) was recorded and used for further calculations.
- f. 0.08 U of PGM and 0.07 U of G6PDH were added in a volume of 1 μl each and the tubes were incubated at 30 °C for 10 min.
- g. The reaction mix of experimental tube was monitored in a spectrophotometer at 340 nm using the reaction mix of reference tube as a reference and the increase in readings noted every minute until constant. The final value of absorbance (*A*₃₄₀, see Calculations, Equations 1, 3) was recorded and used for further calculations.



Vol 6, Iss 18, Sep 20, 2016

2. Soluble starch synthase (SSS) activity assay

Summary of the procedures for SSS and GBSS activity assay is shown in a flowchart (Figure 7).

SSS and GBSS assay

Mix all the components other than ADP-glucose to prepare 900 μ l of SSS or GBSS reaction mixture. Record the volume of enzyme extract added (V_2 for SSS) or the volume of starch granules suspension added (V_7 for GBSS) and the volume of reaction mixture after addition of ADP-glucose (V_4 for both enzymes). Split the volume equally between two Eppendorf tubes. Mark one tube as the reference and another tube as the experiment. Follow section E.2.b of the main protocol for mixture preparation. Each tube should contain 423 μ l of the reaction mixture without ADP-glucose added.




Figure 7. Flowchart of SSS and GBSS assay. Please, refer to the main text of the protocol for specific conditions and buffers/media composition needed for each step of the assay.

www.bio-protocol.org/e1929 Vol 6, Iss 18, Sep 20, 2016

a. For SSS activity assay, two tubes with reaction mixture, each containing 50 mM HEPES-NaOH (pH 7.5), 1.6 mM ADP-glucose*, 1.4 mg amylopectin, 15 mM DTT, and 0.2 ml of enzyme extract in a total volume of 0.45 ml were prepared. The volume of enzyme extract added to starch synthase reaction mixture, as well as total volume of reaction mixture (V_2 and V_4 , correspondingly; see Calculations, Equations 4, 6, 7) were recorded and used for further calculations.

Note*: The reaction mixture is prepared first by mixing all the components other than ADPglucose in a volume of 423 µl as shown in Table 2 (see Recipes, section 3b).

- b. The tubes were pre-incubated in a water bath at 30 °C for 2 min prior to starting the reaction by adding 27 µl of 20 mM ADP-glucose stock solution.
- The reaction in the first tube (reference tube) was stopped immediately after the addition of C. ADP-glucose by incubating the tube in a water bath at 100 °C for 1 min*, whereas the second tube (experimental tube) was incubated at 30 °C for 20 min (t, reaction time, see Calculations, Equations 4, 6, 7) followed by incubating the tube at 100 °C for 1 min*. Subsequently (steps C2d-g), the reactions in the experimental and reference tubes were processed simultaneously in a similar manner.

*Note: Use special microcentrifuge tube lockers (see Materials and Reagents, section D4) to prevent tube opening during boiling. Tube opening at this stage leads to partial evaporation of water and uncontrolled decrease of reaction mixture volume that may seriously affect the accuracy of the measurement.

- d. The reaction mixtures were allowed to cool down to room temperature followed by the centrifuging at 10,000 x g at room temperature for 10 min, 0.3-0.4 ml aliquot of supernatant was collected for further analysis. The volume of the aliquot was recorded and used for further calculations (V_5 , see Calculations, Equations 4, 6, 7).
- 0.2 ml of Solution 1 (see Recipes, section 3c) was added to the aliquot of supernatant from e. step C2d. along with 1.2 U of pyruvate kinase in a volume of 1-2 µl and incubated at 30 °C for 20 min followed by heating the tube at 100 °C for 1 min.
- The reaction mixture was allowed to cool down to room temperature followed by the f. centrifuging at 10,000 x g at room temperature for 10 min, supernatant was transferred to a fresh tube containing 0.4 ml of Solution 2 (see Recipes, section 3d). The final volume of reaction mixture was recorded and used for further calculations (V₃, see Calculations, Equations 4, 6, 7).
- g. 1.4 U hexokinase and 0.35 U G6PDH were added to the final mixture in a volume of 1-2 µl each and the tubes were incubated at 30 °C for 10 min.
- h. The reaction mix of experimental tube was monitored in a spectrophotometer at 340 nm using the reaction mix of reference tube as a reference and the increase in readings noted every minute until constant. The final value of absorbance (A340, see Calculations, Equations 4, 6, 7) was recorded and used for further calculations.

www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

- 3. Granule bound starch synthase (GBSS) activity assay
 - a. Insoluble starch synthase was assayed as described above for SSS (step C2), but the reaction mixture described in section C2a was modified (see Recipes, section 3b). Amylopectin and enzyme extract were omitted, while 0.1 ml aliquot of the starch granules pellet suspension was added as a source of both enzyme and reaction substrate. The suspension was vigorously mixed before being added to the reaction mixture. Settling of granules during the assay was prevented by slow and continuous shaking of the reaction tubes.
 - b. The volume of starch granules pellet suspension added was recorded and used for further calculations (*V*₇, see Calculations, Equation 7).
 - c. All further steps of GBSS assay were similar to that described for SSS (see steps C2b-h).
- 4. Determination of protein content
 - Protein content was determined in grain extracts by Bradford assay (Bradford, 1976) using Bio-Rad Protein Assay Kit and BSA as a standard.
 - b. To perform the analysis, follow the instructions in the manual supplied with the kit available in your laboratory.
 - c. Protein content in the enzyme extract (*C*_{PR}, expressed in *mg ml*⁻¹) was used for calculation of enzymatic activity on the basis of *mg* of protein according to Equations 2, 3, 5, 6 (see Calculations, sections D1c and D2c).
- D. Calculations
 - 1. AGPase
 - a. Enzyme activity of AGPase was expressed in enzyme units. In terms of this protocol, one enzyme unit (*U*) of AGPase is determined as the amount of enzymatic activity that converts 1 micromole of ADP-glucose into glucose-1-phosphate per 1 min at pH 7.5 and at 30 °C.
 - b. Enzyme activity of AGPase on *g* of tissue fresh weight may be calculated according to equation (1).

$$A_{FW} = \frac{Abs_{340}}{6.22 \cdot L} \cdot \frac{V_1 \cdot V_3}{V_2} \cdot \frac{1}{m} \cdot \frac{1}{t} \quad (1)$$

Where,

6.22	-	Micromolar extinction coefficient for NAD(P)H at 340 nm, <i>cm</i> ² µmol ⁻¹ ,	
Abs340	-	Absorbance of reaction mixture of experimental tube measured at 340 nm	
		against the reaction mixture of reference tube,	
AFW	-	Enzyme activity of AGPase calculated on fresh weight basis, $U g^{-1}$,	
L	-	Length of optical path, cm. This parameter is 1.0 cm in most cases,	
т	-	Fresh weight of the sample, g,	
t	-	Reaction time, <i>min</i> ,	

bio-protocol

www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

- V_1 The volume of the extract, *ml*,
- V₂ The volume of the extract added to the reaction mixture, *ml*,
- *V*₃ The volume of final reaction mixture, *ml*.
- c. Enzyme activity of AGPase on *mg* of protein basis may be calculated according to equations (2) and (3).

$$A_{pr} = \frac{A_{FW} \cdot m}{V_1 \cdot C_{pr}} \quad (2)$$

$$A_{pr} = \frac{Abs_{340}}{6.22 \cdot L} \cdot \frac{V_3}{V_2} \cdot \frac{1}{C_{pr}} \cdot \frac{1}{t}$$
(3)

Where,

Apr - Enzyme activity of AGPase calculated on protein basis, U mg⁻¹,

 C_{pr} - Protein concentration in the extract, $mg ml^{-1}$,

and other parameters are the same as for equation (1).

- 2. Starch synthase
 - Enzyme activity of starch synthase (both SSS and GBSS) was expressed in enzyme units.
 In terms of this protocol, one enzyme unit (*U*) of starch synthase (both SSS and GBSS) is determined as the amount of enzymatic activity that converts 1 micromole of ADP-glucose into starch and 1 micromole of ADP per 1 min at pH 7.5 and at 30 °C.
 - Enzyme activity of SSS on g of tissue fresh weight may be calculated according to equation (4).

$$A_{FW} = \frac{Abs_{340}}{6.22 \cdot L} \cdot V_3 \cdot \frac{V_1}{V_2} \cdot \frac{V_4}{V_5} \cdot \frac{1}{m} \cdot \frac{1}{t} \quad (4)$$

Where,

AFW	-	Enzyme activity of SSS calculated on fresh weight basis, $U g^{1}$,
6.22	-	Micromolar extinction coefficient for NAD(P)H at 340 nm, cm ² µmol ⁻¹ ,
Abs340	-	Absorbance of reaction mixture of experimental tube measured at 340 nm
		against the reaction mixture of reference tube,
L	-	Length of optical path, cm. This parameter is 1.0 cm in most cases,
т	-	Fresh weight of the sample, g,
t	-	Reaction time, <i>min,</i>
V_1	-	The volume of the extract, <i>ml</i> ,
V ₂	-	The volume of the extract added to the reaction mixture, <i>ml</i> ,
V ₃	-	The volume of final reaction mixture, <i>ml</i> ,
V4	-	The volume of SSS reaction mixture, <i>ml</i> ,



Vol 6, Iss 18, Sep 20, 2016

- *V*₅ The volume of an aliquot of SSS reaction mixture taken after stopping the reaction by boiling, *ml*.
- c. Enzyme activity of SSS on *mg* of protein basis may be calculated according to equations (5) and (6).

$$A_{pr} = \frac{A_{FW} \cdot m}{V_1 \cdot C_{pr}} \quad (5)$$

$$A_{pr} = \frac{Abs_{340}}{6.22 \cdot L} \cdot \frac{V_3}{V_2} \cdot \frac{V_4}{V_5} \cdot \frac{1}{C_{pr}} \cdot \frac{1}{t}$$
(6)

Where,

Apr - Enzyme activity of SSS calculated on protein basis, U mg⁻¹,

 C_{pr} - Protein concentration in the extract, $mg ml^{-1}$,

and other parameters are the same as for equation (4).

d. Enzyme activity of GBSS on g of tissue fresh weight may be calculated according to equation (7).

$$A_{FW} = \frac{Abs_{340}}{6.22 \cdot L} \cdot V_3 \cdot \frac{V_6}{V_7} \cdot \frac{V_4}{V_5} \cdot \frac{1}{m} \cdot \frac{1}{t}$$
(7)

Where,

AFW	-	Enzyme activity calculated on fresh weight basis,
6.22	-	Micromolar extinction coefficient for NAD(P)H at 340 nm, <i>cm² µmol⁻¹</i> ,
Abs	340 -	Absorbance of reaction mixture of experimental tube measured at 340 nm
		against the reaction mixture of reference tube,
L	-	Length of optical path, cm. This parameter is 1.0 cm in most cases,
т	-	Weight of sample, g,
t	-	Reaction time, <i>min</i> ,
Vз	-	The final volume of the reaction mixture, ml,
V4	-	The volume of GBSS reaction mixture, <i>ml,</i>
V_5	-	The volume of an aliquot of GBSS reaction mixture taken after stopping
		the reaction by boiling, <i>ml,</i>
V_6	-	The volume of starch granules suspension, <i>ml,</i>
V7	-	The volume of the starch granules suspension aliquot added to the
		reaction mixture, <i>ml</i> .
In this	s protoc	ol, we did not measure protein content in starch granule suspension and this

e. In this protocol, we did not measure protein content in starch granule suspension and this is why GBSS activity in the frame of this protocol can be calculated on *g* of fresh weight basis only.

www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

Recipes

- Extraction buffer

 Extraction buffer
 mM HEPES-NaOH (pH 7.5)
 mM MgCl₂
 mM EDTA
 mM DTT
 12.5% (v/v) glycerol
 % (w/v) polyvinylpyrrolidone (PVP-10)
- 2. Stock solutions
 - a. ADP-glucose, disodium salt (MW = 633.31), 20 mM

ADP-glucose stock solution was prepared in the stock container. 0.79 ml of sterile water was added directly to a bottle containing 10 mg of the powder. Solution was aliquoted in small amounts and stored at -20 °C to avoid repetitive thawing.

- b. NADP, sodium salt hydrate (MW = 765.39, anhydrous basis), 10 mM
 7.8 mg of NADP was dissolved in 1 ml sterile water. The solution was aliquoted in small amounts and stored at -20 °C to avoid repetitive thawing. Being frozen, NADP solution is stable for about a year.
- c. PEP monosodium salt hydrate (MW = 190.02, anhydrous basis), 20 mM
 3.9 mg of PEP were dissolved in 1 ml sterile water. The solution was aliquoted in small amounts and stored at -80 °C to avoid repetitive thawing.
- 3. Enzymatic activities reaction mixtures
 - a. AGPase activity reaction mixture (Table 1)

Note: If AGPase activity is low (Abs₃₄₀ is below 0.05), the volume of enzyme extract can be increased at the expense of corresponding decrease of dH₂O volume added.

	Concentration in stock solution, mM	Concentration in reaction mixture, mM	Volume of stock solution per one reaction, µl
HEPES-NaOH, pH 7.5	500	50	45
ADP-glucose	20	1.2	27
Magnesium chloride	500	6.0	5.4
DTT	50	3.0	27
Sodium pyrophosphate	30	5.0	75
dH ₂ O	-	-	220.6
Enzyme extract	-	-	50*
	450		

Table 1. Composition of AGPase activity reaction mixture

www.bio-protocol.org/e1929

Vol 6, Iss 18, Sep 20, 2016

b. Starch synthase activity reaction mixture (Table 2)

Note: If SSS/GBSS activity is low (Abs₃₄₀ is below 0.05), the volume of enzyme extract/starch granule suspension can be increased at the expense of corresponding decrease of dH₂O volume added. In case of SSS, more concentrated DDT can be used. For example, in case of 250 mM stock, 27 μ l instead of 135 μ l of DTT stock has to be added and 308 μ l instead of 200 μ l of enzyme extract can be added to the reaction mixture.

	Concentration in stock solution, mM	Concentration in reaction mixture, mM	Volume of stock one reacti SSS activity assay	-
HEPES-NaOH, pH 7.5	500	50	45	
ADP-glucose	20	1.6	36	
DTT	50	15	135	
Amylopectin (omitted in GBSS assay)	100 mg/ml	3.1 mg/ml	14	-
Enzyme extract (omitted in GBSS assay)	-	-	200*	-
Starch granules suspension (omitted in SSS assay)			-	100*
dH ₂ O	-	-	20	134
Total			450	

Table 2. Composition of starch synthase activity reaction mixture

c. Solution 1 for starch synthase activity assay (Pyruvate kinase reaction mixture) (Table 3) *Note: Solution 1 should be prepared fresh every time while performing the assay.*

Table 3. Composition of Solution 1 (Pyruvate kinase reaction mixture)

	Concentration in stock solution, mM	Concentration in reaction mixture, mM	Volume of stock solution per one reaction, µl
HEPES-NaOH, pH 7.5	500	50	20
Potassium chloride	500	200	80
Magnesium chloride	500	10	4
Phosphoenolpyruvate	20	4	40
dH ₂ O			56
	200		



Vol 6, Iss 18, Sep 20, 2016

d. Solution 2 for starch synthase activity assay (Glucose-6-phosphate dehydrogenase reaction mixture) (Table 4)

Note: Solution 2 should be prepared fresh every time while performing the assay.

	Concentration in stock solution, mM	Concentration in reaction mixture, mM	Volume of stock solution per one reaction, µl
HEPES-NaOH, pH 7.5	500	50	40
Glucose	50	10	80
Magnesium chloride	500	20	16
NADP	10	2	80
dH ₂ O			184
	400		

Table 4. Composition of Solution 2 (G6PDH reaction mixture)

Acknowledgments

This work was supported by a grant from the Natural Sciences and Engineering Research Council of Canada to BTA. The authors would like to acknowledge Ms. Nina Kulichikhina for her technical assistance. This protocol is adapted from Nakamura *et al.* (1989) and Schaffer and Petreikov (1997).

References

- 1. Bradford, M. M. (1976). <u>A rapid and sensitive method for the quantitation of microgram quantities</u> of protein utilizing the principle of protein-dye binding. *Anal Biochem* 72: 248-254.
- Ghosh, H. P. and Preiss, J. (1966). <u>Adenosine diphosphate glucose pyrophosphorylase. A</u> regulatory enzyme in the biosynthesis of starch in spinach leaf chloroplasts. J Biol Chem 241(19): 4491-4504.
- Mukherjee, S., Liu, A., Deol, K. K., Kulichikhin, K., Stasolla, C., Brule-Babel, A. and Ayele, B. T. (2015). <u>Transcriptional coordination and abscisic acid mediated regulation of sucrose transport</u> and sucrose-to-starch metabolism related genes during grain filling in wheat (*Triticum aestivum* <u>L.).</u> *Plant Sci* 240: 143-160.
- 4. Nakamura, Y., Yuki, K., Park, S. Y. and Ohya, T. (1989). <u>Carbohydrate metabolism in the</u> <u>developing endosperm of rice grains</u>. *Plant and Cell Physiology* 30(6): 833-839.
- 5. Schaffer, A. A. and Petreikov, M. (1997). <u>Sucrose-to-starch metabolism in tomato fruit</u> <u>undergoing transient starch accumulation.</u> *Plant Physiology* 113(3): 739-746.
- Vos-Scheperkeuter, G. H., de Boer, W., Visser, R. G., Feenstra, W. J. and Witholt, B. (1986). <u>Identification of granule-bound starch synthase in potato-tubers.</u> *Plant Physiology* 82(2): 411-416.