

Vol 7, Iss 14, Jul 20, 2017 DOI:10.21769/BioProtoc.xxxx

Cell Type-specific Metabolic Labeling of Proteins with Azidonorleucine in Drosophila

Ines Erdmann^{1, 2}, Kathrin Marter^{1, 2, 3}, Oliver Kobler⁴, Sven Niehues^{5, 6}, Julia Bussmann^{5, 6}, Anke Müller^{1, 2}, Julia Abele^{1, 2}, Erik Storkebaum^{5, 6}, Ulrich Thomas^{4, *} and Daniela C. Dieterich^{1, 2, 3, *}

¹Neuronal Plasticity and Communication, Institute for Pharmacology and Toxicology, Otto-von-Guericke-University Magdeburg, Magdeburg, Germany; ²Research Group Neuralomics, Leibniz Institute for Neurobiology, Magdeburg, Germany; ³Center for Behavioral Brain Sciences, Magdeburg, Germany; ⁴Department of Neurochemistry and Molecular Biology, Leibniz Institute for Neurobiology, Magdeburg, Germany; ⁵Molecular Neurogenetics Laboratory, Max Planck Institute for Molecular Biomedicine, Münster, Germany; ⁶Faculty of Medicine, University of Münster, Münster, Germany

*For correspondence: daniela.dieterich@med.ovgu.de; thomas@lin-magdeburg.de

[Abstract] Advanced mass spectrometry technology has pushed proteomic analyses to the forefront of biological and biomedical research. Limitations of proteomic approaches now often remain with sample preparations rather than with the sensitivity of protein detection. However, deciphering proteomes and their context-dependent dynamics in subgroups of tissue-embedded cells still poses a challenge, which we meet with a detailed version of our recently established protocol for cell-selective and temporally controllable metabolic labeling of proteins in *Drosophila*. This method is based on targeted expression of a mutated variant of methionyl-tRNA-synthetase, MetRS^{L262G}, which allows for charging methionine tRNAs with the non-canonical amino acid azidonorleucine (ANL) and, thus, for detectable ANL incorporation into nascent polypeptide chains.

Keywords: Metabolic labeling, Click chemistry, *Drosophila melanogaster*, Proteomic profiling, Protein synthesis

[Background] The protein composition of any given cell is intimately linked to its state of differentiation and functionality. Changes in a cell's proteome may reflect its response to cell-intrinsic cues or to signals originating from elsewhere inside the respective organism or its environment. In turn they are indicative of the significance of those signaling cues. Deciphering proteomes and their dynamics in a cell type-specific fashion has thus become a main focus in current research, reaching a better understanding of molecular events underlying physiological or pathophysiological processes. Any proteomic approach in this direction, however, is challenged by the heterogeneity of cell types that are interconnected within a tissue or organ of interest. In the brain, for instance, different types of neurons and glial cells form the networks required to control animal or human behavior. Moreover, it is well established that information processing within these networks leading to long-term memory is strictly dependent on *de novo* protein synthesis and degradation. While this has been exemplified for a number of neuronal proteins (e.g., immediate early gene proteins), it is obvious that proteins that are



Vol 7, Iss 14, Jul 20, 2017 DOI:10.21769/BioProtoc.xxxx

up- or down-regulated in just a limited number of cells (or even are regulated oppositely in different groups of cells) may easily escape conventional modes of detection, where cellular proteomes are averaged across entire brain areas.

A number of labeling methods for cellular proteomes have been published in the last two decades, e.g., using isotope-coded affinity tags (Gygi et al., 1999) or isobaric tags for relative and absolute quantification (Ross et al., 2004), quantitative proteomic analysis using samples from cells grown in ¹⁴N or ¹⁵N media (Washburn et al., 2002; MacCoss et al., 2003), and stable isotope labeling by amino acids in cell culture (Ong et al., 2002; Andersen et al., 2005). Moreover puromycin (Schmidt et al., 2009) and non-canonical amino acids, e.g., azidohomoalanine (AHA) or homopropargylglycine, in combination with click chemistry have been used to decipher cellular proteomes (Link et al., 2003; Link and Tirrell, 2003; Beatty et al., 2006; Dieterich et al., 2006; Link et al., 2006; Dieterich et al., 2010). All of these strategies, however, fail to uncover cell-type specific proteomes within tissue or organ samples. Most recently, novel strategies to resolve this issue have been reported for C. elegans and Drosophila (Elliott et al., 2014; Erdmann et al., 2015; Yuet et al., 2015). They have in common the use of either a mutated aminoacyl-tRNA synthetase or an orthogonal aminoacyl-tRNA synthetase/tRNA for tagging of newly synthesized proteins with food-supplied non-canonical amino acids. Specifically, we could show that upon cell type-specific expression of a mutant Methionyl-tRNA synthetase (MetRSL262G) as achieved by employing the well-established Gal4/UAS-system, the non-canonical amino acid ANL can be incorporated into proteins of selectable cell types in living Drosophila larvae and adult flies. An accompanying study by Niehues et al. (2015) used this method to show the causal involvement of mutated glycyl-tRNA synthetase in a model for the neurodegenerative Charcot Marie Tooth disease.

ANL-containing proteins can either be analyzed in protein extracts by using biochemistry and mass spectrometry or can be visualized *in situ* by fluorescence microscopy (Erdmann *et al.*, 2015; Niehues *et al.*, 2015). For more information see 'Click Chemistry (CuAAC) and detection of tagged *de novo* synthesized proteins'. The following protocol details the metabolic labeling of proteins in larvae and adult flies with ANL.

Materials and Reagents

- 1. Fly vials (e.g., VWR, catalog number: 734-2254)
- 2. Fly vial plugs (e.g., Carl Roth, catalog number: PK13.1)
- Gal4 activator strains of choice (e.g., C57-Gal4 for muscle-specific expression [from Ulrich Thomas, Magdeburg, Germany], elav^{C155}-Gal4 for pan-neuronal expression [from Bloomington stock center, Bloomington, Indiana, USA], repo-Gal4 for glial expression [from Christian Klämbt, Münster, Germany])
- 4. *UAS-dMetRS*^{L262G} effector strains [available at request from Daniela C. Dieterich & Ulrich Thomas]. As described in Erdmann *et al.* (2015) various lines expressing dMetRS^{L262G} either tagged with 3xmyc or EGFP are available



Vol 7, Iss 14, Jul 20, 2017 DOI:10.21769/BioProtoc.xxxx

Note: We traditionally use ONM. The standard corn meal medium has also been successfully used in Niehues et al. (2015) for ANL labeling, thus, we anticipate that other media can be used as well without any limitations.

- 5. Otto-normal-medium (ONM, see Recipes)
 - a. Agar-Agar (Carl Roth, catalog number: 5210)
 - b. Semolina (local food store)
 - c. Mashed raisins (local food store)
 - d. Baker's yeast (local food store)
 - e. Sugar beet sirup (local food store)
 - f. Honey (local food store)
 - g. Tap water
 - h. 20% (w/v) Nipagin (see Recipes)
 - i. Methyl-4-hydroxybenzoate (Merck, catalog number: 106757)
 - ii. Propyl-4-hydroxybenzoate (Merck, catalog number: 107427)
 - iii. 100 % ethanol (Th. Geyer, catalog number: 2246)
- 6. 200 mM ANL stock solution (for the synthesis of ANL see [Link et al., 2007; Ngo et al., 2009; Erdmann et al., 2015]) (see Recipes)

Equipment

- 1. Beaker (kitchen/household grade)
- 2. Immersion blender (kitchen/household grade)
- 3. Paintbrush (art supplies)
- 4. Fly incubator (e.g., SANYO, model: MIR-553)
- 5. Hotplate (kitchen/household grade)
- 6. Pot (kitchen/household grade)
- 7. Tablespoon (kitchen/household grade)

Procedure

- A. Preparation of ANL-containing fly food medium (Figure 1)
 - 1. Thaw baker's yeast.
 - 2. Add semolina and Agar-Agar to 0.33 L water. Stir from time to time until swelling is completed.
 - 3. Heat mashed raisins, yeast, sugar beet syrup and honey in 0.66 L water. Boil the mixture for 5 min while stirring constantly.
 - 4. Add the semolina-Agar-Agar-mixture and boil once again. Don't forget to stir constantly, as the mixture might braise at the bottom of the pot.
 - 5. Cool down the ONM to 50 °C. Stir every 15-20 min for 1 min.



Vol 7, Iss 14, Jul 20, 2017 DOI:10.21769/BioProtoc.xxxx

- 6. Add the Nipagin and stir for at least 1 min until the Nipagin is homogenously distributed in the food.
- 7. Add 2 ml of ANL stock solution to 100 ml ONM in a beaker for a final concentration of 4 mM ANL. Mix for 1 min using an immersion blender.
- 8. Aliquot ANL-containing ONM (2-4 ml) into fly vials and let cool down completely at room temperature.
- 9. Plug the vials.
- 10. ANL-containing ONM can be stored for approximately two weeks at 4 °C. Discard vials once the ONM detaches from the vial wall.



Figure 1. Preparation of ANL-containing fly food medium. A. Ingredients for ONM are shown: R (mashed raisins); S (semolina); BS (sugar beet syrup); H (honey). B. Semolina and Agar-Agar are added to water and allowed to swell, SAA. C. Raisins and Baker's yeast (Y), sugar beet syrup (BS), and honey (H), are added into water and boiled for 5 min. D. Afterwards, the SAA mixture is added and boiled once more. After cooling down to 50 °C with stirring from time to time, Nigapin and ANL are added and mixed thoroughly. E and F. Media is then aliquoted and allowed to cool down completely before storage at 4 °C.

- B. Cell type-specific expression of dMetRS^{L262G}-variants in *Drosophila* larvae and flies
 - 1. Collect virgin female flies of the respective activator strain.



Vol 7, Iss 14, Jul 20, 2017 DOI:10.21769/BioProtoc.xxxx

Note: The number of flies and hence the number of offspring depends on your experimental design including the type of analysis (e.g., MS/MS) and the effectiveness and degree of cell selectivity of the Gal4-activator in use.

2. Cross virgin female flies of the activator strain to male flies of a UAS-dMetRS^{L262G}-effector strain (Erdmann *et al.*, 2015).

Note: Crosses can as well be set up reciprocally, that is, you may use virgin female flies of a UAS-dMetRS^{L262G}-effector strain and cross them to males of your activator strain.

- C. Cell type-specific labeling of proteins with ANL-examples
 - 1. Long-term labeling of proteins in larval body wall muscles:
 - a. Raise appropriate crosses on ANL-containing ONM such that ANL is present during all developmental stages of the progeny.
 - b. Mid- to late 3rd larval stage is reached after approximately 5 days when crosses are raised at 25 °C and after approximately 10-12 days when kept at 18 °C.
 - c. Perform dissection of larval body walls according to 'Click Chemistry (CuAAC) and detection of tagged de novo synthesized proteins'. See also Bellen and Budnik (2000). Note: Biochemical approaches on whole larval extracts have proven difficult to perform time and again, perhaps due to lytic activities. We therefore recommend separating body walls (mainly comprising muscles, epithelia, cuticle, trachea and sensory neurons) from all other tissues.
 - 2. Long-term labeling of proteins in fly heads (e.g., if dMetRS^{L262G} is expressed in the CNS, compound eyes and/or antenna):
 - a. Raise crosses on ANL-containing ONM so that ANL is present throughout development of the progeny.
 - b. Remove parental flies before eclosure of the offspring.
 - c. Prepare fly heads from adult progeny according to 'Click Chemistry (CuAAC) and detection of tagged *de novo* synthesized proteins'.
 - 3. Short-term labeling of proteins in larval body wall muscles or brains:
 - a. Raise crosses on ANL-free ONM for 1-2 days at 25 °C.
 - b. Transfer parental flies onto fresh ANL-free ONM for 4-6 h. Remove parental flies.
 - c. Keep vials for 72 ± 2 h at 25 °C.
 - d. Wash early 3rd instar larvae out of the food with warm tap water and rinse them into a mesh basket.
 - e. Transfer 3rd instar larvae onto ANL-containing ONM using a paintbrush.
 - f. Keep larvae at 25 °C for 24 h.
 - g. Prepare the larval body walls or larval brains according to 'Click Chemistry (CuAAC) and detection of tagged *de novo* synthesized proteins'.



Vol 7, Iss 14, Jul 20, 2017 DOI:10.21769/BioProtoc.xxxx

- 4. Short-term labeling of proteins in fly heads:
 - Raise appropriate crosses on ANL-free ONM either at 25 °C or 18 °C depending on your experimental design.
 - b. Discard parental flies before progeny ecloses.
 - c. Transfer progeny flies on ANL-containing ONM for a period of time according to your experimental design. For certain analyses it may also be considered to place flies back onto ANL-free food.
 - d. Prepare fly heads according to 'Click Chemistry (CuAAC) and detection of tagged *de novo* synthesized proteins'.

Data analysis

ANL incorporation into fly or larval proteins can be analyzed after performing copper-catalyzed azide-alkyne cycloaddition (CuAAC, 'click chemistry') as described in Erdmann *et al.* (2015) and the accompanying bio-protocol. General fly and larvae viabilities upon ANL incorporation can be analyzed by assessing *e.g.*, locomotor behavior and/or hatching rates as described in Erdmann *et al.* (2015).

Notes

- Reproducibility and variability: Uptake of ANL and thus subsequent labeling efficiency may vary
 when heterogeneous animal numbers are raised in the vials, as the food will be differently
 mashed through depending on the number of larvae. Also of note, larvae will show in general
 stronger ANL incorporation compared to adult flies.
- 2. ANL-containing ONM is used best within two weeks after preparation to yield reproducible results.

Recipes

- 1. 20% (w/v) Nipagin
 - 150 g methyl-4-hydroxybenzoate
 - 50 g propyl-4-hydroxybenzoate
 - Dissolve in 100% (abs) ethanol
- 2. Otto-normal-medium (ONM)
 - 8.3 g Agar-Agar
 - 50 g semolina
 - 40 g mashed raisins
 - 2 cubes of yeast (42 g each)
 - 1 tablespoon sugar beet syrup



Vol 7, Iss 14, Jul 20, 2017 DOI:10.21769/BioProtoc.xxxx

1 tablespoon honey 6.6 ml 20% (w/v) Nipagin

200 mM ANL stock solution
 34.4 mg per 1 ml ddH₂O

Acknowledgments

This work was supported by the German Research Foundation (DFG) with an SFB 779 grant to D.C.D., SFB 854 grants to U.T. and D.C.D. and a Leibniz Society PAKT grant (LGS SynaptoGenetics) to D.C.D. and U.T. This protocol was adapted from Erdmann *et al.*, 2015.

References

- 1. Andersen, J. S., Lam, Y. W., Leung, A. K., Ong, S. E., Lyon, C. E., Lamond, A. I. and Mann, M. (2005). Nucleolar proteome dynamics. *Nature* 433(7021): 77-83.
- Beatty, K. E., Liu, J. C., Xie, F., Dieterich, D. C., Schuman, E. M., Wang, Q. and Tirrell, D. A. (2006). <u>Fluorescence visualization of newly synthesized proteins in mammalian cells.</u> Angew Chem Int Ed Engl 45(44): 7364-7367.
- 3. Bellen, H. J. and Budnik, V. (2000). *Drosophila*, a laboratory manual. In: Ashburner, M., Hawley, S. and Sullivan. B (Eds). Cold Spring Harbor Laboratory. chap. 11.
- 4. Dieterich, D. C., Hodas, J. J., Gouzer, G., Shadrin, I. Y., Ngo, J. T., Triller, A., Tirrell, D. A. and Schuman, E. M. (2010). *In situ* visualization and dynamics of newly synthesized proteins in rat hippocampal neurons. *Nat Neurosci* 13(7): 897-905.
- 5. Dieterich, D. C., Link, A. J., Graumann, J., Tirrell, D. A. and Schuman, E. M. (2006). <u>Selective identification of newly synthesized proteins in mammalian cells using bioorthogonal noncanonical amino acid tagging (BONCAT)</u>. *Proc Natl Acad Sci U S A* 103(25): 9482-9487.
- Elliott, T. S., Townsley, F. M., Bianco, A., Ernst, R. J., Sachdeva, A., Elsasser, S. J., Davis, L., Lang, K., Pisa, R., Greiss, S., Lilley, K. S. and Chin, J. W. (2014). <u>Proteome labeling and protein identification in specific tissues and at specific developmental stages in an animal.</u> *Nat Biotechnol* 32(5): 465-472.
- 7. Erdmann, I., Marter, K., Kobler, O., Niehues, S., Abele, J., Muller, A., Bussmann, J., Storkebaum, E., Ziv, T., Thomas, U. and Dieterich, D. C. (2015). <u>Cell-selective labelling of proteomes in *Drosophila melanogaster*. *Nat Commun* 6: 7521.</u>
- 8. Gygi, S. P., Rist, B., Gerber, S. A., Turecek, F., Gelb, M. H. and Aebersold, R. (1999). Quantitative analysis of complex protein mixtures using isotope-coded affinity tags. *Nat Biotechnol* 17(10): 994-999.
- 9. Link, A. J., Mock, M. L. and Tirrell, D. A. (2003). Non-canonical amino acids in protein engineering. *Curr Opin Biotechnol* 14(6): 603-609.



Vol 7, Iss 14, Jul 20, 2017 DOI:10.21769/BioProtoc.xxxx

- 10. Link, A. J. and Tirrell, D. A. (2003). <u>Cell surface labeling of Escherichia coli via copper(I)-catalyzed [3+2] cycloaddition</u>. *J Am Chem Soc* 125(37): 11164-11165.
- Link, A. J., Vink, M. K., Agard, N. J., Prescher, J. A., Bertozzi, C. R. and Tirrell, D. A. (2006).
 <u>Discovery of aminoacyl-tRNA synthetase activity through cell-surface display of noncanonical amino acids.</u> *Proc Natl Acad Sci U S A* 103(27): 10180-10185.
- 12. Link, A. J., Vink, M. K. and Tirrell, D. A. (2007). <u>Synthesis of the functionalizable methionine surrogate azidohomoalanine using Boc-homoserine as precursor.</u> *Nat Protoc* 2(8): 1884-1887.
- 13. MacCoss, M. J., Wu, C. C., Liu, H., Sadygov, R. and Yates, J. R., 3rd (2003). <u>A correlation algorithm for the automated quantitative analysis of shotgun proteomics data.</u> *Anal Chem* 75(24): 6912-6921.
- Ngo, J. T., Champion, J. A., Mahdavi, A., Tanrikulu, I. C., Beatty, K. E., Connor, R. E., Yoo, T. H., Dieterich, D. C., Schuman, E. M. and Tirrell, D. A. (2009). <u>Cell-selective metabolic labeling of proteins</u>. *Nat Chem Biol* 5(10): 715-717.
- Ong, S. E., Blagoev, B., Kratchmarova, I., Kristensen, D. B., Steen, H., Pandey, A. and Mann, M. (2002). <u>Stable isotope labeling by amino acids in cell culture</u>, <u>SILAC</u>, as a simple and accurate approach to expression proteomics. *Mol Cell Proteomics* 1(5): 376-386.
- Ross, P. L., Huang, Y. N., Marchese, J. N., Williamson, B., Parker, K., Hattan, S., Khainovski, N., Pillai, S., Dey, S., Daniels, S., Purkayastha, S., Juhasz, P., Martin, S., Bartlet-Jones, M., He, F., Jacobson, A. and Pappin, D. J. (2004). Multiplexed protein quantitation in Saccharomyces cerevisiae using amine-reactive isobaric tagging reagents. Mol Cell Proteomics 3(12): 1154-1169.
- 18. Schmidt, E. K., Clavarino, G., Ceppi, M. and Pierre, P. (2009). <u>SUnSET, a nonradioactive</u> method to monitor protein synthesis. *Nat Methods* 6(4): 275-277.
- Washburn, M. P., Ulaszek, R., Deciu, C., Schieltz, D. M. and Yates, J. R., 3rd (2002). <u>Analysis of quantitative proteomic data generated via multidimensional protein identification technology</u>. *Anal Chem* 74(7): 1650-1657.
- 20. Yuet, K. P., Doma, M. K., Ngo, J. T., Sweredoski, M. J., Graham, R. L., Moradian, A., Hess, S., Schuman, E. M., Sternberg, P. W. and Tirrell, D. A. (2015). <u>Cell-specific proteomic analysis in Caenorhabditis elegans.</u> *Proc Natl Acad Sci U S A* 112(9): 2705-2710.