

Isolation and Long-term Cultivation of Mouse Alveolar Macrophages

Clara Jana-Lui Busch^{1, 2}, Jérémy Favret^{1, 2, 3}, Laufey Geirsdóttir³, Kaaweh Molawi^{2, 3}
and Michael H. Sieweke^{1, 2, 3, *}

¹Center for Regenerative Therapies Dresden (CRTD), Technische Universität Dresden, Fetscherstraße 105, Dresden, Germany; ²Max-Delbrück-Centrum für Molekulare Medizin in der Helmholtzgemeinschaft (MDC), Berlin, Germany; ³Aix Marseille University, CNRS, INSERM, CIML, Marseille, France

*For correspondence: Michael.sieweke@tu-dresden.de

[Abstract] Alveolar macrophages (AM) are tissue-resident macrophages that colonize the lung around birth and can self-maintain long-term in an adult organism without contribution of monocytes. AM are located in the pulmonary alveoli and can be harvested by washing the lungs using the method of bronchoalveolar lavage (BAL). Here, we compared different conditions of BAL to obtain high yields of murine AM for *in vitro* culture and expansion of AM. In addition, we describe specific culture conditions, under which AM proliferate long-term in liquid culture in the presence of granulocyte-macrophage colony-stimulating factor. This method can be used to obtain large numbers of AM for *in vivo* transplantation or for *in vitro* experiments with primary mouse macrophages.

Keywords: Macrophage, Alveolar, Lungs, Self-renewal, Bronchoalveolar lavage, Primary cell culture

[Background] AM are resident tissue macrophages of the lungs with critical importance for immune regulation and surfactant homeostasis (Kopf *et al.*, 2015). Due to their localization in the airspace of the alveoli, AM are directly exposed to inhaled air and pathogens or other aerosolized particles. Consequently, AM play a crucial role in the initiation or suppression of inflammatory responses and are the subject of investigation in numerous studies that explore mechanisms of pulmonary diseases (Hodge *et al.*, 2007; Sun and Metzger, 2008; Happle *et al.*, 2014; Schneider *et al.*, 2014a; Machiels *et al.*, 2017; Yu *et al.*, 2017). Interestingly, AM in the mouse originate from fetal monocytes and are able to self-maintain their numbers *in vivo* without contribution of bone-marrow-derived monocytes under steady state conditions (Guilliams *et al.*, 2013; Hashimoto *et al.*, 2013). AM are unique in that they reside outside of the body surface and are directly exposed to the external environment. They can therefore be isolated with minimal tissue disturbance using bronchoalveolar lavage (BAL). We previously demonstrated that the self-renewal property of AM harvested by BAL is maintained in culture by growing AM long-term in liquid media or serially re-plating AM in semi-solid media (Soucie *et al.*, 2016; Imperatore *et al.*, 2017). Here, we describe the methodological advancements in BAL and specific culturing conditions optimized for long-term culture and high AM yields.

To obtain sufficient quality and numbers of AM, we tested different cell harvesting conditions and developed a culture method that allows long-term maintenance of AM *in vitro*. Our BAL method was based on earlier studies using pre-warmed PBS and EDTA for detaching AM from lung alveoli (Steele *et al.*, 2003; Zhang *et al.*, 2008). Additionally, serum was added for cellular protection during the isolation

period. Importantly, both omitting EDTA in the lavage buffer and using pre-cooled PBS for lavage resulted in lower AM yield in our hands. Our comparisons demonstrated that the variable with the largest effect on yield was the temperature of the BAL buffer (Figure 1). Whereas many earlier studies obtain BAL cell numbers below 2×10^5 per wild-type mouse (Table 1), our comparative analysis indicated that BAL cell numbers can be considerably increased by these optimizations. We generally obtained up to 5×10^5 - 7×10^5 live (Trypan-Blue-negative) cells per mouse using this method (Figure 1). FACS analysis revealed that more than 98% of cells were alive at the time of recording using the Zombie Violet Fixable Viability Kit, independent of whether the BAL was performed with 4 °C PBS or 37 °C PBS/EDTA/FBS. Our experience has shown that when performing dozens of BALs on the same day long waiting times on ice will result in higher cell death unless low amounts of FBS are added to the BAL buffer (typically 0.5%, although we have good experience with up to 2%).




Figure 1. Comparison of BAL conditions using pre-cooled or pre-warmed PBS with or without 2 mM EDTA. Either 4 °C PBS without EDTA, 4 °C PBS with 2 mM EDTA and 0.5% FBS, or 37 °C PBS with 2 mM EDTA and 0.5% FBS was used. Numbers show the total amount of living cells (Trypan-Blue-negative) per BAL treatment per mouse. Each symbol denotes the mean cell count of 3 technical replicates of an individual mouse; horizontal lines indicate the mean, error bars show standard error of the mean (SEM); one-way ANOVA with Tukey's multiple comparisons test; ns, non-significant.

Table 1. Comparison of the efficiency of different BAL protocols. Total cell number in BAL fluid obtained in cited studies. Only counts from WT animals (or comparable conditions, such as control-treated WT animals) were considered. WT denotes wild-type mouse.

Mouse strain	Total cell number	Volume	Buffer conditions	References
WT (4 weeks)	20,000	2 x 0.4 ml	0.5 mM EDTA	Yu <i>et al.</i> , 2017
WT (control)	40,000	2 x 0.8 ml		Yuan <i>et al.</i> , 2015
WT (floxed)	30-40,000	3 x 0.4 ml		Schneider <i>et al.</i> , 2014b
WT (control)	50,000	2 x 1 ml	EDTA, 4 °C	Machiels <i>et al.</i> , 2017
WT	60,000	3 x 0.4 ml		Schneider <i>et al.</i> , 2014a
WT (6 weeks)	50-150,000	3 x 1 ml	0.5 mM EDTA	van de Laar <i>et al.</i> , 2016
WT (1 year)	100-200,000	3 x 1 ml	0.5 mM EDTA	van de Laar <i>et al.</i> , 2016
WT (control)	250,000	5 x 1 ml	5 mM EDTA	Qian <i>et al.</i> , 2018
WT (1 year)	250-300,000	5 x 1 ml		Suzuki <i>et al.</i> , 2014
WT	300-500,000	3 x 2 ml	0.5 mM EDTA, 37 °C	Zhang <i>et al.</i> , 2008
WT	400-500,000	3 x 1 ml	2% BSA, 2 mM EDTA	Dong <i>et al.</i> , 2018
WT (6-8 weeks)	500-700,000	10 x 1 ml	0.5% FBS, 2 mM EDTA, 37 °C	Soucie <i>et al.</i> , 2016, Imperatore <i>et al.</i> , 2017 and this study

Increased efficiency in harvesting BAL cells is advantageous for *in vivo* reconstitution of multiple recipient animals such as strains devoid of endogenous AM, such as GM-CSF-receptor-deficient mice (Guilliams *et al.*, 2013; van de Laar *et al.*, 2016). High starting cell numbers and viability of harvested cells will accelerate establishment of long-term AM cultures and improve cellular yield. We also noticed several culture conditions outlined in the detailed protocol that affect the quality, yield, doubling time and durability of the cultures. In order to avoid cell activation or death, several parameters need to be controlled. Firstly, we used exclusively sterile supplies and applied sterile handling techniques to avoid activation of AM. Secondly, we use non-treated plastic ware (not tissue-culture treated plates or dishes) and a gentle detachment protocol. Together, these technical improvements will be helpful for starting and maintaining a long-term AM culture.

Biochemical and genetic manipulations of macrophages that require a large number of cells could so far only be done in cell lines, such as RAW 264.7 or J774A.1 cells (Ralph and Ninkovic, 1975; Raschke *et al.*, 1978), oncogene-transformed cells, for example Myc (Baumbach *et al.*, 1986; Baumbach *et al.*, 1987) or SV40 transformed macrophages like IC-21 (Walker and Demus, 1975) or MH-S cells (Mbawuke and Herscovitz, 1989), macrophages differentiated from progenitors (Zhang *et al.*, 2008; Fejer *et al.*, 2013) or non-transformed but genetically modified macrophages, such as Maf-DKO macrophages (Aziz *et al.*, 2009). The ability to obtain large numbers of normal unmodified AM in culture allows such experiments in primary resident macrophages. Our studies on macrophage self-renewal mechanisms serve as an example (Soucie *et al.*, 2016; Imperatore *et al.*, 2017).

Materials and Reagents

1. 15-ml conical tubes (Corning, catalog number: 352196)

2. Bottle-top vacuum filter with 0.22 µm membrane (Corning, catalog number: 431161)
3. Plastic storage bottle (Corning, catalog number: 430281)
4. 70-µm sterile cell strainer (BD, catalog number: 340633)
5. 1-ml syringe (Braun, catalog number: 9161406V)
6. 18-G cannula (Braun, catalog number: 4667123)
7. Petri dish 94/16 mm (Greiner Bio-one, catalog number: 633181)
8. Non-treated 6-well plate (Nunc™, catalog number: 150239)
9. C57BL/6 mice (aged 6-10 weeks)
10. PBS, pH 7.2 (Thermo Fisher Scientific, Gibco™, catalog number: 20012019)
11. EDTA stock solution (*e.g.*, 0.5 M, pH 8.0)
12. Hemolysis buffer (self-made or commercial, *e.g.*, Morphisto, catalog number: 12146)
13. Trypan Blue solution 0.4% (Sigma-Aldrich, catalog number: T8154)
14. RPMI 1640 Medium, no glutamine (Thermo Fisher Scientific, Gibco™, catalog number: 31870025)
15. Fetal bovine serum (Testing of different batches is recommended)
16. Gentamicin sulphate 50 mg/ml in aqueous solution (Lonza, catalog number: BE02-012E)
17. Penicillin-Streptomycin (10,000 U/ml) (Thermo Fisher Scientific, Gibco™, catalog number: 15140122)
18. Sodium Pyruvate (100 mM) (Thermo Fisher Scientific, Gibco™, catalog number: 11360070)
19. GlutaMAX™ Supplement (Thermo Fisher Scientific, Gibco™, catalog number: 35050038)
20. Conditioned medium from J558L cell line transfected with murine GM-CSF cDNA as a source for GM-CSF (Zal *et al.*, 1994; Stockinger *et al.*, 1996; Rayasam, 2015)
21. ESGRO Complete Accutase (Merck, catalog number: SF006)
22. EGTA stock solution (*e.g.*, 0.5 M, pH 8.0)
23. UltraComp eBeads™ Compensation Beads (Thermo Fisher Scientific, Invitrogen, catalog number: 01-2222-41)
24. Zombie Violet Fixable Viability Kit (BioLegend, catalog number: 423113)
25. FACS antibodies (as indicated in Table 2)
26. BAL buffer (see Recipes)
27. Complete medium (see Recipes)
28. AM culture medium (see Recipes)
29. Detachment medium (see Recipes)

Equipment

1. Pipettes
2. Mouse dissection tools (scissors, forceps)
3. Water bath set to 37 °C
4. Refrigerated benchtop centrifuge for spinning conical tubes

5. Hemocytometer (Roth, catalog number: T729.1)
6. Incubator (37 °C, 5% CO₂)
7. Inverse microscope

Procedure

A. Harvest alveolar macrophages by bronchoalveolar lavage (BAL)

1. For each mouse, prepare a 15-ml conical tube filled with 3 ml complete medium (see Recipes).
2. Warm-up BAL buffer (see Recipes) to 37 °C in a water bath. Keep warm during the whole procedure.
3. Euthanize the mouse by cervical dislocation without rupturing the jugular vein or the trachea to avoid exposing AM to CO₂ or isoflurane, which could affect functional properties of AM.
4. Using dissection tools, remove the skin, ribcage and muscles to expose both lungs and trachea. Avoid cutting or rupturing blood vessels.

Note: Since methods for the surgical exposure of lungs and the trachea have been published previously in this journal, the reader is referred to those protocols for instructions (Han and Ziegler, 2013; Tibbitt and Coquet, 2016; Jhingran et al., 2016; Sun et al., 2017).

5. Use a fine scissor to make a small incision in the upper part of the trachea just below the larynx. The part of the trachea facing downwards (away from the experimenter) should remain intact, do not cut through the whole trachea.
6. Use the incision to insert a slightly blunted 18-G cannula and direct the cannula 5 mm deeper into the trachea towards the lungs. Take care not to damage lung tissue.
7. Attach a 1 ml syringe filled with 1 ml warm BAL buffer onto the inserted cannula.
8. Inject 1 ml buffer while fixating the cannula position with the other hand.
9. Pull the plunger to collect BAL fluid in the syringe. About 800-900 µl can be recovered. Observe that the pressure should not be too high, otherwise the alveoli will burst and BAL fluid will be lost. Upon injection and collection, the lungs should visibly inflate and deflate.
10. Filter collected BAL fluid through a 70 µm cell strainer into the 15-ml tube with 3 ml complete medium from Step A1.
11. Repeat Steps A6-A10 for 9 more times each time with fresh warm BAL buffer. Pool cells into the same 15-ml tube.
12. Collect cells by centrifugation at 300 x g, 5 min at 4 °C. Remove supernatant. The cell pellet should be white. A red/pink color indicates that blood was accidentally collected during the BAL.
13. Add 1 ml hemolysis buffer for 2 min incubation at room temperature to lyse residual erythrocytes. Fill up tube with complete medium to stop lysis and collect cells by centrifugation as before. Remove supernatant. The color of the cell pellet should be white now.
14. Resuspend cell pellet in 500 µl BAL buffer and take a sample for counting using a hemocytometer chamber after staining with Trypan Blue to exclude dead cells. Count only live (Trypan-Blue negative) cells.

15. Calculate the total cell number per BAL. Typically, 5×10^5 - 7×10^5 live cells per adult wild-type mouse aged 6-8 weeks can be recovered when using pre-warmed BAL buffer, containing PBS with 2 mM EDTA and 0.5% serum.
16. Proceed to cell staining and flow cytometry analysis or *in vitro* cultivation.

B. Flow cytometric analysis of alveolar macrophages

1. Block unspecific binding sites on cells with TruStain fcX and concomitantly stain with Zombie Violet in 200 μ l cold PBS (without FBS) at 4 °C in the dark for 15 min (see Table 2).
2. Wash cells with cold BAL buffer by centrifugation at 300 x *g* for 5 min at 4 °C.
3. Stain cells in a volume of 100 μ l per 1 million cells according to Table 2 for 30 min at 4 °C in the dark using BAL buffer.

Table 2. FACS reagents used for staining BAL AM

	Conjugate	Clone	Dilution	Company, Cat. No.
TruStain fcX	–	93	1:200	BioLegend, 101320
ZombieViolet	BV421	–	1:400	BioLegend, 423113
SiglecF	PE	E50-2440	1:200	BD, 552126
CD11c	PE/Cy7	N418	1:400	BioLegend, 117318

4. Prepare compensation beads for each antibody conjugate.
5. Wash cells with BAL buffer, resuspend in 200 μ l BAL buffer for recording.
6. Record cells by flow cytometry after acquiring the compensation beads. AM are double-positive for SiglecF and CD11c (Figures 2A-2C), and > 98% viable (Figure 2D).




Figure 2. FACS analysis of BAL AM. A. Simple gating strategy for exclusion of doublets, dead cells. AM are SiglecF⁺ and CD11c⁺ positive. B-C. BAL cells harvested with pre-warmed BAL buffer containing EDTA are phenotypically not different from BAL cells harvested using pre-cooled PBS. Each symbol denotes 1 mouse. Typically, > 95% of BAL cells are AM. D. Viability analysis of BAL singlet cells assessed by staining with Zombie Violet Fixable Dye.

C. Cultivation of alveolar macrophages

1. Collect cells by centrifugation as before. Remove supernatant.
2. Plate 3×10^5 - 4×10^5 cells per well of a non-treated 6-well plate in 3 ml pre-warmed AM culture medium (see Recipes).

Note: Typically, 3×10^5 - 4×10^5 cells are plated in 1 well of a 6-well plate. If BAL cells of several mice are pooled, 1.1-1.2 million cells can be plated in a non-treated 94 mm Petri dish in 10 ml pre-warmed AM culture medium.

3. Add gentamicin to the AM culture (1:1,000).

Note: Gentamicin is omitted after the first medium change.

4. Incubate at 37 °C, 5% CO₂.
5. Replace culture supernatant after 6-18 h with fresh AM culture medium and discard the supernatant.

Note: Cells will adhere fully within a few hours after the first plating and we do not keep cells in suspension at the first medium exchange. However, for subsequent medium exchanges, the cells in the supernatant are collected as well since a typical AM culture consists of both adherent and suspended cells (see also Notes section).

6. Change medium every 2 days until the cell culture reaches confluency.
7. To change medium, transfer the medium and suspension cells into a 15 ml-tube. Add 2 ml warm AM culture medium to the well with adherent cells to prevent drying-out. Collect the suspension

cells using centrifugation at 300 x g 5 min. Resuspend the pelleted cells in 1 ml warm AM culture medium and combine with adherent cells.

Note: Freshly harvested primary AM will double every 7-10 days (Soucie et al., 2016). If the majority of AM appear stretched (spindle-like) and activated, increasing the amount of conditioned medium or adding recombinant GM-CSF might help; however, proliferative capacity will be limited and it might advisable to start a new culture (see Figure 3 for exemplary images of early AM culture).




Figure 3. Representative images of AM culture within the first days after plating the cells.

A. AM culture with predominantly round-shaped cells that are partly floating and partly adherent on Day 1 after plating. B. Same culture as (A) on Day 2. C. Same culture as (A) on Day 4. D. Example for an AM culture with a large fraction of elongated, dark cells on Day 4 after plating. Arrowheads indicate dividing cells, 100× magnification.

8. To detach cells from a confluent well, collect suspension cells into a 15-ml tube.
9. Add 750 µl detachment medium (see Recipes) to 1 well of a 6-well plate (or 3 ml to a 94 mm Petri dish) and incubate for 10-30 min at 37 °C.

Note: AM are very adherent and prone to rupture when using too harsh detachment procedures. Thus, the use of non-treated plastic ware and proper detachment medium is important (see Recipes). Ruptured cells in the culture medium might affect both activation status and proliferative capacity of AM.

10. Detachment of cells can be supported by pipetting on the plastic bottom gently to avoid cellular damage (see Notes).
11. Pool detached cells with cells in suspension and centrifuge cells at 300 x g for 5 min.

12. Resuspend cell pellet in 1 ml warm AM cultured medium and take a sample for counting using a hemocytometer chamber after staining with Trypan Blue to exclude dead cells. Count only live (Trypan-Blue negative) cells.
13. If the cell number has doubled, add 5 ml warm AM culture medium and split into 2 wells of a 6-well plate (or correspondingly to 2x 94 mm Petri dishes). In general, the cell number plated is maintained around the values indicated above in the note to Step C2.

Note: Earlier, we could show that AM culture remains proliferative for at least 10 passages (Soucie et al., 2016). Since then, we have experience with AM cultures that remain proliferative even beyond 20 passages with no indication of a decline in proliferative capacity.

Data analysis

Harvested cells were counted manually using a hemocytometer and considering only Trypan-Blue-negative cells. Stained cells were recorded on a BD LSRFortessa with 5 lasers using BD FACSDiva software and analyzed using FlowJo v10. Microscope images were acquired on an inverse microscope (Leica DMI1) equipped with a digital camera (MC120). Gating was performed as indicated in Figure 2A. To test for statistically significant differences between the means of three groups (Figure 1), one-way ANOVA with Tukey's multiple comparisons test was performed using GraphPad Prism 7. No data points were excluded.

Notes

1. Proliferative AM are round-shaped and semi-adherent. Re-plating of suspension AM will result in part of the cells attaching to the new well, while the other part remains in suspension. Take care to not lose the suspended cells when changing medium as this will reduce the number of proliferative cells and slow the expansion of the culture.
2. When detaching cells, do not pipet the cell suspension up and down extensively, this might affect the viability of the culture; if cells do not detach readily, collect detachment medium containing already detached cells and perform another round of incubation with fresh detachment medium and/or increase the incubation time. Late-passage cells require shorter incubation times (~ 5 min) than early-passage cultures (up to 30 min).
3. The percentage of conditioned medium should be titrated after preparation of each batch by testing the growth of AM in the presence of various amounts of conditioned medium (e.g., 1, 2, 5, 10% in complete medium). In our batches, we use typically 1%-3% conditioned medium diluted into complete medium (e.g., 100 µl in 10 ml), which corresponds roughly to 2-5 ng/ml purified recombinant mouse GM-CSF.
4. We successfully replaced conditioned medium with 20 ng/ml recombinant GM-CSF (Peprotech) for long-term culture. Lower GM-CSF concentrations might be sufficient but have not been tested.

Recipes

1. BAL buffer
PBS
2 mM EDTA (dilute 1:250 from 0.5 M EDTA stock solution)
0.5% Fetal bovine serum (FBS)
Sterile-filter using vacuum filtration and keep at 4 °C until use
2. Complete medium
RPMI 1640
1x GlutaMAX
1x Pyruvate
1x Penicillin/Streptomycin
10% FBS
Sterile-filter using vacuum filtration and keep at 4 °C until use
3. AM culture medium
Supplement complete medium with 1-5% conditioned medium containing mouse GM-CSF (needs to be titrated)
Pre-warm an aliquot in a water bath to 37 °C before use
4. Detachment medium
ESGRO Complete Accutase
1 mM EGTA (dilute 1:500 from 0.5 M EGTA stock solution)
Aliquot and freeze at -20 °C. Pre-warm an aliquot in the water bath to 37 °C before use

Acknowledgments

We thank Stephanie Vargas Aguilar and Sethuraman Subramanian for technical discussions and helpful comments on the manuscript, Philippe Pierre for J558L cells, and the flow cytometry and animal caretaking facilities of the Max-Delbrück-Centrum for assistance. This study was supported by institutional grants from Institut National de la Santé et de la Recherche Médicale, Centre National de la Recherche Scientifique, and Aix-Marseille University to the CIML and grants to M.H. Sieweke from the Agence Nationale de la Recherche (ANR-11-BSV3-0026 and ANR-17-CE15-0007-01), the Institut National du Cancer (InCA13-10/405/AB-LC-HS), Fondation pour la Recherche Médicale (DEQ. 20110421320), and the European Research Council (ERC) under the European Union's Horizon 2020 research and innovation program (grant agreement number 695093 MacAge). K. Molawi was supported by a Human Frontier Science Program long-term fellowship and Stiftung Charité, L. Geiersdottir by the Fondation ARC pour la Recherche sur le Cancer. M.H. Sieweke is a Berlin Institute of Health–Einstein visiting fellow at MDC and an Alexander von Humboldt Professor at TU Dresden.

Competing interests

The authors declare no competing financial interests.

Ethics

Animal husbandry and mouse work were conducted in accordance with the German Animal Welfare legislation, after the approval by the Landesamt für Gesundheit und Soziales (for work in Berlin, following the guidelines of the Institutional Animal Care and Use Committee of the Max Delbrück Centrum für Molekulare Medizin) and after the approval by the Landesdirektion Sachsen (for work in Dresden).

References

1. Aziz, A., Soucie, E., Sarrazin, S. and Sieweke, M. H. (2009). [MafB/c-Maf deficiency enables self-renewal of differentiated functional macrophages](#). *Science* 326(5954): 867-871.
2. Baumbach, W. R., Keath, E. J. and Cole, M. D. (1986). [A mouse c-myc retrovirus transforms established fibroblast lines *in vitro* and induces monocyte-macrophage tumors *in vivo*](#). *J Virol* 59(2): 276-283.
3. Baumbach, W. R., Stanley, E. R. and Cole, M. D. (1987). [Induction of clonal monocyte-macrophage tumors *in vivo* by a mouse c-myc retrovirus: rearrangement of the CSF-1 gene as a secondary transforming event](#). *Mol Cell Biol* 7(2): 664-671.
4. Dong, Y., Poon, G. F. T., Arif, A. A., Lee-Sayer, S. S. M., Dosanjh, M. and Johnson, P. (2018). [The survival of fetal and bone marrow monocyte-derived alveolar macrophages is promoted by CD44 and its interaction with hyaluronan](#). *Mucosal Immunol* 11(3): 601-614.
5. Fejer, G., Wegner, M. D., Gyory, I., Cohen, I., Engelhard, P., Voronov, E., Manke, T., Ruzsics, Z., Dolken, L., Prazeres da Costa, O., Branzk, N., Huber, M., Prasse, A., Schneider, R., Apte, R. N., Galanos, C. and Freudenberg, M. A. (2013). [Nontransformed, GM-CSF-dependent macrophage lines are a unique model to study tissue macrophage functions](#). *Proc Natl Acad Sci U S A* 110(24): E2191-2198.
6. Guilliams, M., De Kleer, I., Henri, S., Post, S., Vanhoutte, L., De Prijck, S., Deswarte, K., Malissen, B., Hammad, H. and Lambrecht, B. N. (2013). [Alveolar macrophages develop from fetal monocytes that differentiate into long-lived cells in the first week of life via GM-CSF](#). *J Exp Med* 210(10): 1977-1992.
7. Han, H. and Ziegler, S. F. (2013). [Bronchoalveolar lavage and lung tissue digestion](#). *Bio-protocol* 3(16): e859.
8. Happel, C., Lachmann, N., Skuljec, J., Wetzke, M., Ackermann, M., Brenning, S., Mucci, A., Jirmo, A. C., Groos, S., Mirenska, A., Hennig, C., Rodt, T., Bankstahl, J. P., Schwerk, N., Moritz, T. and Hansen, G. (2014). [Pulmonary transplantation of macrophage progenitors as effective and](#)

- [long-lasting therapy for hereditary pulmonary alveolar proteinosis. *Sci Transl Med* 6\(250\): 250ra113.](#)
9. Hashimoto, D., Chow, A., Noizat, C., Teo, P., Beasley, M. B., Leboeuf, M., Becker, C. D., See, P., Price, J., Lucas, D., Greter, M., Mortha, A., Boyer, S. W., Forsberg, E. C., Tanaka, M., van Rooijen, N., Garcia-Sastre, A., Stanley, E. R., Ginhoux, F., Frenette, P. S. and Merad, M. (2013). [Tissue-resident macrophages self-maintain locally throughout adult life with minimal contribution from circulating monocytes. *Immunity* 38\(4\): 792-804.](#)
10. Hodge, S., Hodge, G., Ahern, J., Jersmann, H., Holmes, M. and Reynolds, P. N. (2007). [Smoking alters alveolar macrophage recognition and phagocytic ability: implications in chronic obstructive pulmonary disease. *Am J Respir Cell Mol Biol* 37\(6\): 748-755.](#)
11. Imperatore, F., Maurizio, J., Vargas Aguilar, S., Busch, C. J., Favret, J., Kowenz-Leutz, E., Cathou, W., Gentek, R., Perrin, P., Leutz, A., Berruyer, C. and Sieweke, M. H. (2017). [SIRT1 regulates macrophage self-renewal. *EMBO J* 36\(16\): 2353-2372.](#)
12. Jhingran, A., Kasahara, S. and Hohl, T. M. (2016). [Flow cytometry of lung and bronchoalveolar lavage fluid cells from mice challenged with fluorescent *Aspergillus* reporter \(FLARE\) conidia. *Bio-protocol* 6\(18\): e1927.](#)
13. Kopf, M., Schneider, C. and Nobs, S. P. (2015). [The development and function of lung-resident macrophages and dendritic cells. *Nat Immunol* 16\(1\): 36-44.](#)
14. Machiels, B., Dourcy, M., Xiao, X., Javaux, J., Mesnil, C., Sabatel, C., Desmecht, D., Lallemand, F., Martinive, P., Hammad, H., Guillailliams, M., Dewals, B., Vanderplasschen, A., Lambrecht, B. N., Bureau, F. and Gillet, L. (2017). [A gammaherpesvirus provides protection against allergic asthma by inducing the replacement of resident alveolar macrophages with regulatory monocytes. *Nat Immunol* 18\(12\): 1310-1320.](#)
15. Mbawuike, I. N. and Herscowitz, H. B. (1989). [MH-S, a murine alveolar macrophage cell line: morphological, cytochemical, and functional characteristics. *J Leukoc Biol* 46\(2\): 119-127.](#)
16. Qian, G., Jiang, W., Zou, B., Feng, J., Cheng, X., Gu, J., Chu, T., Niu, C., He, R., Chu, Y. and Lu, M. (2018). [LPS inactivation by a host lipase allows lung epithelial cell sensitization for allergic asthma. *J Exp Med* 215\(9\): 2397-2412.](#)
17. Ralph, P. and Nakoinz, I. (1975). [Phagocytosis and cytolysis by a macrophage tumour and its cloned cell line. *Nature* 257\(5525\): 393-394.](#)
18. Raschke, W. C., Baird, S., Ralph, P. and Nakoinz, I. (1978). [Functional macrophage cell lines transformed by Abelson leukemia virus. *Cell* 15\(1\): 261-267.](#)
19. Rayasam, A. (2015). [Isolating and culturing dendritic cells \(Dendritic Cell J558 Protocol\). *Protoc Exch* 22-24.](#)
20. Schneider, C., Nobs, S. P., Heer, A. K., Kurrer, M., Klinke, G., van Rooijen, N., Vogel, J. and Kopf, M. (2014a). [Alveolar macrophages are essential for protection from respiratory failure and associated morbidity following influenza virus infection. *PLoS Pathog* 10\(4\): e1004053.](#)
21. Schneider, C., Nobs, S. P., Kurrer, M., Rehrauer, H., Thiele, C. and Kopf, M. (2014b). [Induction of the nuclear receptor PPAR- \$\gamma\$ by the cytokine GM-CSF is critical for the differentiation of fetal](#)

- [monocytes into alveolar macrophages](#). *Nat Immunol* 15(11): 1026-1037.
22. Soucie, E. L., Weng, Z., Geirsdottir, L., Molawi, K., Maurizio, J., Fenouil, R., Mossadegh-Keller, N., Gimenez, G., VanHille, L., Beniazza, M., Favret, J., Berruyer, C., Perrin, P., Hacohen, N., Andrau, J. C., Ferrier, P., Dubreuil, P., Sidow, A. and Sieweke, M. H. (2016). [Lineage-specific enhancers activate self-renewal genes in macrophages and embryonic stem cells](#). *Science* 351(6274): aad5510.
23. Steele, C., Marrero, L., Swain, S., Harmsen, A. G., Zheng, M., Brown, G. D., Gordon, S., Shellito, J. E. and Kolls, J. K. (2003). [Alveolar macrophage-mediated killing of *Pneumocystis carinii* f. sp. muris involves molecular recognition by the Dectin-1 \$\beta\$ -glucan receptor](#). *J Exp Med* 198(11): 1677-1688.
24. Stockinger, B., Zal, T., Zal, A. and Gray, D. (1996). [B cells solicit their own help from T cells](#). *J Exp Med* 183: 891-899.
25. Sun, F., Xiao, G. and Qu, Z. (2017). [Murine bronchoalveolar lavage](#). *Bio-protocol* 7(10): e2287.
26. Sun, K. and Metzger, D. W. (2008). [Inhibition of pulmonary antibacterial defense by interferon- \$\gamma\$ during recovery from influenza infection](#). *Nat Med* 14(5): 558-564.
27. Suzuki, T., Arumugam, P., Sakagami, T., Lachmann, N., Chalk, C., Salles, A., Abe, S., Trapnell, C., Carey, B., Moritz, T., Malik, P., Lutzko, C., Wood, R. E. and Trapnell, B. C. (2014). [Pulmonary macrophage transplantation therapy](#). *Nature* 514(7523): 450-454.
28. Tibbitt, C. and Coquet, J. M. (2016). [House dust mite extract and cytokine instillation of mouse airways and subsequent cellular analysis](#). *Bio-protocol* 6(14): e1875.
29. van de Laar, L., Saelens, W., De Priek, S., Martens, L., Scott, C. L., Van Isterdael, G., Hoffmann, E., Beyaert, R., Saeys, Y., Lambrecht, B. N. and Guillems, M. (2016). [Yolk sac macrophages, fetal liver, and adult monocytes can colonize an empty niche and develop into functional tissue-resident macrophages](#). *Immunity* 44(4): 755-768.
30. Walker, W. S. and Demus, A. (1975). [Antibody-dependent cytolysis of chicken erythrocytes by an *in vitro*-established line of mouse peritoneal macrophages](#). *J Immunol* 114(2 pt 2): 765-769.
31. Yu, X., Buttgerit, A., Lelios, I., Utz, S. G., Cansever, D., Becher, B. and Greter, M. (2017). [The cytokine TGF- \$\beta\$ promotes the development and homeostasis of alveolar macrophages](#). *Immunity* 47: 903-912.e4.
32. Yuan, X., Shan, M., You, R., Frazier, M. V., Hong, M. J., Wetsel, R. A., Drouin, S., Seryshev, A., Song, L. Z., Cornwell, L., Rossen, R. D., Corry, D. B. and Kheradmand, F. (2015). [Activation of C3a receptor is required in cigarette smoke-mediated emphysema](#). *Mucosal Immunol* 8(4): 874-885.
33. Zal, T., Volkmann, A. and Stockinger, B. (1994). [Mechanisms of tolerance induction in major histocompatibility complex class II-restricted T cells specific for a blood-borne self-antigen](#). *J Exp Med* 180(6): 2089-2099.
34. Zhang, X., Goncalves, R. and Mosser, D. M. (2008). [The isolation and characterization of murine macrophages](#). *Curr Protoc Immunol* Chapter 14: Unit 14.1.