

Analysis of *Pseudomonas aeruginosa* c-di-GMP High and Low Subpopulations Using Flow-assisted Cell Sorting (FACS) and Quantitative Reverse Transcriptase PCR (qRT-PCR)

Catherine R. Armbruster¹ and Matthew R. Parsek^{2, *}


¹Department of Microbiology and Molecular Genetics, University of Pittsburgh, Pittsburgh, PA, USA;

²Department of Microbiology, University of Washington, Seattle, WA, USA

*For correspondence: parsem@uw.edu

[Abstract] Cyclic diguanylate monophosphate (c-di-GMP) is a second messenger signaling molecule that drives the transition from planktonic to the biofilm mode of growth in many bacterial species. *Pseudomonas aeruginosa* has at least two surface sensing systems that produce c-di-GMP in response to surface attachment, the Wsp and Pil-Chp systems. We recently used a plasmid-based c-di-GMP reporter (pP_{cdrA}::gfp) to describe how the Wsp system generates heterogeneity in surface sensing, resulting in two physiologically distinct subpopulations of cells during early biofilm formation. One subpopulation has elevated c-di-GMP and produces biofilm matrix, serving as the founders of initial microcolonies. The other subpopulation has low c-di-GMP and engages in surface motility, allowing for exploration of the surface. Here, we describe the protocol for a key experiment to confirm our initial observation of c-di-GMP heterogeneity during surface sensing: the use of flow-assisted cell sorting (FACS) to isolate subpopulations of cells with high and low c-di-GMP reporter activity, followed by quantitative Reverse Transcriptase PCR (qRT-PCR) of genes that are known to be transcriptionally regulated in response to cellular c-di-GMP levels (*pelA*, *pslA*). This protocol can be adapted by others to isolate subpopulations of high- and low- c-di-GMP *P. aeruginosa* cells that are genetically identical, but phenotypically distinct for future experiments examining specific mRNA transcripts as we did or, presumably, for additional applications like RNAseq, proteomics, or TNseq.

Graphical abstract:



Keywords: c-di-GMP, Flow cytometry, Single-cell, c-di-GMP reporter, Biofilms

[Background] The second messenger signaling molecule cyclic diguanylate monophosphate (c-di-GMP) allows bacteria to rapidly modify their cell surface in response to environmental conditions.

Diguanylate cyclases (DGC) are proteins containing a GGDEF amino acid sequence motif that synthesize c-di-GMP from two molecules of GTP, and phosphodiesterases (PDE) are proteins with EAL or HD-GYP motifs that hydrolyze c-di-GMP into linear pGpG or GMP (Hengge, 2009). The genomes of *Pseudomonas aeruginosa* lab strains PAO1 and PA14 each encode almost 40 of these putative c-di-GMP-modulating enzymes, the majority of which contain additional predicted domains involved in environmental sensing (including PAS, GAF, and HAMP domains) and signal transduction (e.g., REC domains) (Kulasakara *et al.*, 2006). Through the activity of these DGCs and PDEs that modulate intracellular levels of c-di-GMP, *P. aeruginosa* integrates diverse environmental signals in order to coordinate a variety of bacterial behaviors, including the transition from planktonic to the biofilm mode of growth.

In *P. aeruginosa*, elevated c-di-GMP levels promote cell surface adhesiveness by interacting with specific receptors and effectors that activate production of biofilm extracellular matrix components and suppress expression of motility genes. The c-di-GMP binding transcriptional regulator FleQ acts in complex with FleN to inversely regulate transcription of genes involved in flagellar biosynthesis and biofilm matrix production (Arora *et al.*, 1997; Hickman and Harwood, 2008; Baraquet *et al.*, 2012; Baraquet and Harwood 2016). When cellular levels of c-di-GMP are high, such as during growth on a surface, FleQ binds c-di-GMP and derepresses transcription of genes encoding the Pel and Psl polysaccharide biosynthetic machinery, as well as the adhesin CdrA (Borlee *et al.*, 2010). The c-di-GMP reporters described in this protocol contain the promoter region of *cdrA* and can be used in *P. aeruginosa* strains that express functional FleQ and FleN proteins (Figure 1).

Phenotypic heterogeneity is a common phenomenon in nature, occurring when a population of genetically identical organisms (in the present case, bacteria) display diverse behaviors, despite experiencing the same environmental conditions (Ackermann, 2015). Testing hypotheses related to changes in gene expression in phenotypically heterogeneous populations requires the ability to separate out these sub-populations. For example, when analyzed as a single batch (combining both sub-populations together), the high gene expression by one sub-population would be averaged with the low expression from the other sub-population, resulting in the erroneous conclusion that expression of the gene of interest is at some intermediate level. Therefore, the purpose of developing this protocol was to separate genetically identical, but apparently physiologically distinct sub-populations of *P. aeruginosa* cells with high and low intracellular levels of c-di-GMP in order to determine whether these two subpopulations may differ in their expression of genes involved in biofilm formation (genes encoding Pel and Psl polysaccharide biosynthetic machineries). The protocol uses either of two versions of a plasmid-based, transcriptional reporter of c-di-GMP levels in *P. aeruginosa* (pP_{cdrA}::gfp and pP_{cdrA}::gfp_{ASV}) (Armbruster *et al.*, 2019). These two reporter plasmids differ by which GFP allele they express, either the stable GFP (GFPmut3) or an unstable variant with a shorter half-life (GFP_{ASV}) (Andersen *et al.*, 1998; Rybtke *et al.*, 2012). The benefit of the stable GFP version of the reporter is that it tends to have a higher dynamic range, whereas the benefit of the short half-life GFP version is that it can be used in time course experiments when the levels of c-di-GMP may decrease. The protocol also makes use of a number of controls to help draw flow cytometry gates that separate cells with high and low c-di-GMP reporter

activity (as measured by GFP fluorescence). Controls for gating include promotorless vector control plasmids harboring stable or unstable GFP (pMH487 or pMH489) and 2 *P. aeruginosa* strains with clean deletions of key proteins in the Wsp chemosensory system (PAO1 $\Delta wspR$ and PAO1 $\Delta wspF \Delta pelC \Delta pslBCD$) that will have constitutively low or high levels of c-di-GMP reporter activity (Güvener and Harwood, 2007).




Figure 1. The plasmid-based transcriptional reporters, $pP_{cdrA}::gfp$ and $pP_{cdrA}::gfp_{ASV}$, are sensitive to the intracellular concentration of c-di-GMP in *P. aeruginosa*. Under conditions of low intracellular c-di-GMP (top), FleQ and FleN repress the promoter for the adhesin, CdrA (*PcdrA*), and no GFP is produced. Under conditions of high c-di-GMP (bottom), such as when *P. aeruginosa* is attached to a surface, c-di-GMP binds to FleQ, the *cdrA* promoter is de-repressed, and GFP is made. Yellow star = c-di-GMP. RBSII is an enhanced ribosomal binding site, T0 and T1 are transcriptional terminators, and *Gent^R* and *Amp^R* represent resistance genes for gentamicin and ampicillin, respectively. This is the Seattle c-di-GMP reporter that was originally published in Rybtke *et al.* (2012).

The Wsp system is a “surface sensing” system in *P. aeruginosa* with homology to bacterial chemosensory signal transduction complexes, that produces c-di-GMP in response to surface contact (Hickman *et al.*, 2005). The Wsp system senses a yet to be identified signal related to surface contact through WspA, a membrane-bound receptor homologous to methyl-accepting chemotaxis proteins (MCPs), transducing a signal to the CheA homolog, WspE. Upon autophosphorylation, the sensor kinase WspE activates the diguanylate cyclase WspR through phosphorylation of a conserved aspartate residue on WspR’s receiver domain. Phosphorylation of WspR leads to the formation of large subcellular WspR-P aggregates that synthesize c-di-GMP (Huangyutitham *et al.*, 2013). WspC (a methyltransferase) and WspF (a methylesterase) are involved in adaptation of the Wsp system to stimulation. The PAO1

$\Delta wspR$ deletion strain is incapable of producing c-di-GMP through the Wsp system and is a helpful control strain for gating a “low c-di-GMP” population. In contrast, the PAO1 $\Delta wspF \Delta peIC \Delta psIBCD$ strain has constitutively high c-di-GMP produced through the Wsp system due to the deletion of the *wspF* methylesterase gene. This strain also has deletions in genes required for Pel and Psl polysaccharide biosynthesis ($\Delta peIC \Delta psIBCD$), to abrogate clumping due to polysaccharide overproduction in the PAO1 $\Delta wspF$ background.

Materials and Reagents

1. 1 ml serological pipette
2. 2 ml tube
3. Control bacterial strains (to determine low GFP and high GFP gates)
4. Experimental bacterial strains

This protocol utilizes *P. aeruginosa* strains harboring the following plasmids that are available upon request from the laboratory of Dr. Matthew R. Parsek: pP_{cdrA}::*gfp*, pP_{cdrA}::*gfp*_{ASV}, pMH487, pMH489. However, in principle, the protocol can be adapted for use with other fluorescent reporters.

5. Genomic DNA isolated from *P. aeruginosa* PAO1
6. Gentamicin
7. LB agar
8. Paraformaldehyde (PFA; Millipore Sigma, catalog number: 158127)
9. TRIzol LS (Thermo Fisher Scientific, Waltham, MA)
10. RQ1 DNase I (Promega, Madison, WI)
11. iTaq Universal SYBR Green One-Step kit (Bio-Rad, Hercules, CA)
12. LB broth (see Recipes)
 - a. Peptone (Fisher Scientific, catalog number: BP1420)
 - b. Gibco™ Bacto™ Yeast extract (Fisher Scientific, catalog number: DF0127-17-9)
 - c. Sodium chloride (NaCl; Fisher Scientific, catalog number: S271-3)
13. Modified Jensen's Chemically Defined Media (see Recipes)
 - a. Sodium chloride (NaCl; Fisher Scientific, catalog number: S271-3)
 - b. Potassium phosphate dibasic anhydrous (K₂HPO₄; Fisher Scientific, catalog number: P288)
 - c. L-Glutamic acid monosodium salt hydrate (Millipore Sigma, catalog number: G5889)
 - d. L-Valine (Fisher Scientific, catalog number: AC140811000)
 - e. L-Phenylalanine (Fisher Scientific, catalog number: AC130311000)
 - f. Dextrose (D-Glucose; Fisher Scientific, catalog number: D16)
 - g. Magnesium sulfate heptahydrate (MgSO₄; Fisher Scientific, catalog number: M63)
 - h. Calcium chloride dihydrate (CaCl₂; Fisher Scientific, catalog number: C79)
 - i. Iron (II) sulfate heptahydrate (FeSO₄; Millipore Sigma, catalog number: F7002)
 - j. Zinc sulfate heptahydrate (ZnSO₄; Millipore Sigma, catalog number: Z4750)

14. Sterile PBS (see Recipes)

- Sodium chloride (NaCl; Fisher Scientific, catalog number: S271-3)
- Potassium chloride (KCl; Fisher Scientific, catalog number: BP366-500)
- Sodium phosphate dibasic (Na₂HPO₄; Fisher Scientific, catalog number: S375)
- Potassium phosphate monobasic (KH₂PO₄; Fisher Scientific, catalog number: P285)

Equipment

- 80 °C freezer
- Micropipettes
- Benchtop microcentrifuge
- FACSAriaII (BD Biosciences, San Jose, CA) or equivalent flow cytometer that is capable of measuring forward and side scattered light and GFP fluorescence, as well as cell sorting
- CFX96 Touch Real-Time PCR detection system (Bio-Rad, Hercules, CA) or equivalent
- NanoDrop or equivalent spectrophotometer that is capable of measuring optical density at the following wavelengths: 230, 260, and 280 nanometers (nm)

Software

- FlowJo (FlowJo LLC, Ashland, OR) or equivalent software to analyze flow cytometry datasets

Procedure

The following procedure can be adapted for use with cells grown under a variety of c-di-GMP inducing conditions prior to FACS analysis. The methods described below were previously used in a study to examine heterogeneity in c-di-GMP production during surface sensing and early biofilm formation by *P. aeruginosa* (Armbruster *et al.*, 2019).

- Prepare control and experimental bacterial strains for inoculation into the experimental condition (Figure 2)
 - Streak *P. aeruginosa* reporter strains from a -80 °C glycerol freezer stock onto LB agar with 100 µg/ml of gentamicin. Incubate plates overnight at 37 °C.
 - Prepare an overnight liquid culture by resuspending a single colony into 3 ml LB broth with 100 µg/ml of gentamicin. Grow the liquid culture shaking, overnight, at 37 °C.
 - Generate an inoculum of *P. aeruginosa* cells grown to mid-log in planktonic culture. Back dilute overnight, stationary phase cells 1:1000 in a final volume of 3 ml LB with 100 µg/ml gentamicin. Continue to incubate the back-diluted culture at 37 °C with shaking for 2 h. The purpose of generating a mid-log planktonic culture of cells is to start the experiment under conditions where c-di-GMP levels (and GFP fluorescence) should be lowest in the cells.

4. After 2 h, begin checking the optical density of your culture at 600 nanometers (OD₆₀₀) until cultures have reached mid-log phase (OD₆₀₀ ~0.5).




Figure 2. Graphical representation of Steps in Part A. Inocula of reporter and control strains are prepared to mid-logarithmic phase prior to inoculating into the experimental conditions in Part B.

- B. Grow reporter and control strains under c-di-GMP inducing experimental condition. Set up an additional culture growing to mid-log as a control for single cells
 1. Centrifuge 3 ml of mid-log cells for 1 min at $> 6,800 \times g$ at room temperature (15-25 °C).
 2. Resuspend concentrated cells in a final volume of 1 ml LB with 100 µg/ml gentamicin, to achieve a final OD₆₀₀ of 0.8.
 3. Incubate resuspended cells statically in a sterile, glass culture tube for 4 h at 37 °C. The glass culture tube acts as the substratum for *P. aeruginosa* to attach to and initiate c-di-GMP signaling over 4 h.
 4. Set up the single-cell control by further back diluting one of the strains from Step A4 to OD₆₀₀ 0.05 in LB and growing shaking at 37 °C for 4 h.
- C. Harvest surface-attached cells and fix with paraformaldehyde (PFA)
 1. Prepare the single-cell control sample from Step B4 by centrifuging 1 ml of mid-log cells for 1 min at $> 6,800 \times g$ at room temperature (15-25 °C). Rinse twice with PBS, then resuspend the cell pellet in 1 ml of ice-cold 6% PFA. Store this control tube on wet ice.
 2. Use a 1 ml serological pipette to gently remove LB containing the planktonic fraction of each experimental culture from Step B3.

3. Rinse the inside of the tube with 2 ml of ice cold PBS. Perform 3 rinses. Store tubes on wet ice between washes.
4. Use a 1 ml serological pipette to transfer 1 ml of ice cold, 6% PFA into each tube. Use the pipette to scrape the sides of the tube and vortex thoroughly to resuspend surface-attached cells.
5. Incubate resuspended samples on ice for 20 min for PFA fixation.
Note: We used PFA to fix the bacterial cells instead of RNAlater because RNAlater is known to reduce the fluorescence intensity of GFP (Zaitoun et al., 2010).
6. Within an hour of harvesting samples and immediately prior to FACS analysis, centrifuge to pellet PFA-fixed cells, wash once in ice cold PBS, and resuspend in 1 ml ice cold PBS. Temporarily store samples on wet ice until they are used for FACS.

D. Analyze populations by flow cytometry and draw gates for sorting

We used a FACSAriaII (BD Biosciences, San Jose, CA) flow cytometer with a 70 μ m nozzle to sort *P. aeruginosa* cells harboring the c-di-GMP reporter into high and low c-di-GMP populations for downstream qRT-PCR. The FACSAriaII was cleaned with 10% bleach followed by PBS for 90 min before and after all runs.

1. Use the following controls to draw gates for high and low GFP levels (Figure 3):
 - a. Single cell control: Wild type PAO1 cells grown planktonically to mid-log. Use this sample first to draw a gate around single cells, using forward and side scatter. This will allow you to avoid counting cell debris or large aggregates of cells as single events when sorting.
 - b. GFP fluorescence controls: Use these to draw gates where you expect to see high and low GFP fluorescence from the *PcdrA* reporter (*i.e.*, cells with high or low c-di-GMP levels, respectively).
 - i. If using the regular half-life GFP reporter ($pP_{cdrA}::gfp$):
 - 1) Constitutively high GFP controls:
PAO1 $\Delta wspF \Delta peiC \Delta ps/BCD$ $pP_{cdrA}::gfp$
PAO1 Tn7::P(A1/04/03)::GFPmut
 - 2) Low GFP control:
PAO1 $\Delta wspR$ $pP_{cdrA}::gfp$
 - 3) Constitutively “off” reporter controls:
PAO1 pMH487
Wild type PAO1
 - ii. If using the short half-life GFP reporter ($pP_{cdrA}::gfp^{ASV}$):
 - 1) Constitutively high GFP control:
PAO1 $\Delta wspF \Delta peiC \Delta ps/BCD$ $pP_{cdrA}::gfp^{ASV}$
 - 2) Low GFP control:
PAO1 $\Delta wspR$ $pP_{cdrA}::gfp^{ASV}$
 - 3) Constitutively “off” reporter control:
PAO1 pMH489




Figure 3. Controls for setting up GFP gates for flow cytometry with the c-di-GMP reporter.

Brackets indicate gates for 'on' (GFP above 1.7×10^2 RFU) or 'off' (GFP below 1.7×10^2 RFU) reporter cells and the number above the bracket indicates the percentage of cells that fall within that gate. (A) Wild type PAO1 cells were used to determine the background level of fluorescence on the BD Aria III for GFP measurements. This population of cells falls below 10^3 RFU. (B) A *P. aeruginosa* strain constitutively expressing stable GFP (PAO1 Tn7::P(A1/04/03)::GFPmut) was used to determine gating for cells with the highest amount of GFP, with the population ranging from 10^3 to 10^5 RFU. (C) Surface grown PAO1 $\Delta wspF \Delta pel \Delta psI$ harboring the $pP_{cdrA}::gfp_{ASV}$ was used to validate the gate for collection of cells with high reporter activity (10^3 to 10^5 RFU; 91.6% of the population). (D) Example of gating for reporter 'on' cells from wild type PAO1 $pP_{cdrA}::gfp_{ASV}$ cells that had been attached to glass in LB medium for 4 h. Approximately 56.6% of the population falls into the reporter 'on' population. (E) Example of gating for 'on' cells from wild type PA14 $pP_{cdrA}::gfp_{ASV}$ cells that had been attached to glass in LB medium for 4 h prior to FACS sorting. Approximately 12.8% of the population falls into the reporter 'on' population. (F) Example of gating for 'on' cells from wild type PA14 $pP_{cdrA}::gfp_{ASV}$ cells that had been attached to glass in Jensen's media (a condition in which Pel is more abundantly produced than in LB) for 4 h. Approximately 31.9% of the population falls into the reporter 'on' population. This figure was reproduced from Figure 1-figure supplement 5 in Armbruster *et al.* (2019).

E. Sort cells directly into TRIzol LS

1. Prepare tubes to sort cells into by adding 750 μ l of Trizol LS for every 250 μ l of volume that you plan to sort into that tube (*i.e.*, sorted bacterial cells suspended in PBS sheath fluid). Pre-chill tubes with Trizol LS at 4 $^{\circ}$ C.

2. Sort cells at 4 °C by flow assisted cell sorting (FACS) to collect 100,000 events (single cells) into TRIzol LS. Monitor the final volume to maintain a 3:1 ratio of Trizol LS to sorted cell volume and add additional Trizol LS as needed.
3. After sorting is completed on each tube, vortex the samples for 10 s and store on wet ice.

F. Extract RNA and test for gDNA contamination

1. Extract RNA immediately following FACS-sorting.
 - a. Incubate samples for 4 min at 90 °C. Invert to mix.
 - b. Transfer samples to ice.
 - c. Add 200 µl of ice cold chloroform and shake vigorously.
 - d. Centrifuge at 6,800 x g for 15 min at 4 °C.
 - e. Remove 500 µl of the colorless aqueous layer and transfer to a clean 2 ml tube.
 - f. Add the following to the tube:
 - 1 µl glycogen
 - 40 µl of 3 M NaOAc, pH 5.0
 - 1,200 µl of 100% ethanol
 - g. Vortex to mix and store at -20 °C for 2 h or overnight.
 - h. Centrifuge RNA samples for 30 min at 14,000 x g, 4 °C.
 - i. Remove supernatant and wash RNA pellet with 800 µl 70% ethanol.
 - j. Centrifuge for 30 min, remove the supernatant, and briefly air dry the pellet at room temperature.
 - k. Resuspend the pellet in 20 µl RNase-free water.
 - l. Heat at 65 °C for 5 min and vortex to fully resuspend the pellet.
 - m. Proceed with RQ1 DNase I digestion.
2. Digest any contaminating DNA present in the RNA samples with RQ1 DNase I (Promega, Madison, WI) following the manufacturer's protocol.
3. Check for genomic DNA contamination by PCR to detect the *rplU* gene. Include PAO1 genomic DNA as a positive control.
 - a. Primers (amplifies a 307 bp product):
 - rplU*-F: 5'-CGCAGTGATTGTTACCGGTG-3'
 - rplU*-R: 5'-AGGCCTGAATGCCGGTGATC-3'
 - b. Set up PCR reactions using your preferred Polymerase following manufacturer's instructions (Table 1):

Table 1. Example PCR recipe using Phusion polymerase

	Per reaction (μl)
5x HF buffer	4
100% DMSO	0.6
10mM dNTPs	0.4
100 μM <i>rplU</i> -F	0.04
100 μM <i>rplU</i> -R	0.04
ddH ₂ O	14.68
RNA sample	0.5
Phusion polymerase	0.2
Total volume:	20

- c. Program your thermocycler appropriately to allow for amplification of an approximately 300 bp amplicon with the polymerase you are using (Table 2).

Table 2. Example PCR conditions to amplify *rplU* when using Phusion polymerase

Temperature (°C)	Time
98	30 s
98	10 s
54	15 s
72	20 s
72	5 min
12	∞

Repeat 34 times

- d. Run amplicons on an electrophoresis gel to check for gDNA contamination. Include a DNA ladder and the positive control sample. Contaminated RNA samples will have a band around 300 bp in length, the same length as the positive control (PAO1 gDNA)'s amplicon. For any samples with gDNA contamination, repeat the treatment with Promega RQ1 DNase.
- G. Assess expression of target genes by quantitative Reverse Transcriptase PCR (qRT-PCR)
- The following steps describe how to measure expression of 2 target genes (*peIA* and *psIA*) and one control gene (*ampR*) by quantitative Reverse Transcriptase PCR (qRT-PCR) using the iTaq Universal SYBR Green One-Step kit (Bio-Rad, Hercules, CA) and a CFX96 Touch Real-Time PCR detection system (Bio-Rad, Hercules, CA). Conditions will vary depending on the target genes, qPCR kit, and the thermalcycler used.

1. Determine total RNA concentration and purity. Use a NanoDrop or other spectrophotometer to measure optical density of each RNA sample at 260 nm, 230 nm, and 280 nm. Dilute samples 1:50 in 100 µl H₂O and repeat the measurement.
 - a. To estimate concentration, an OD₂₆₀ of 1 unit equals approximately 40 µg/ml RNA.
 - b. A sufficiently pure RNA sample will have 260/280 and 260/230 ratios above 1.8.
2. Prepare genomic DNA standards from 20ng to 0.02pg using PAO1 gDNA diluted in TE buffer.
3. Set up qRT-PCR reactions using the manufacturer's instructions. Use 2 or 3 µl of RNA per reaction, depending on the final concentration of RNA you measured above (add no more than 500 ng of RNA to each reaction).
 - a. Control primers:

ampR-F-qPCR: 5'-GCGCCATCCCTTCATCG-3'

ampR-R-qPCR: 5'-GATGTGCGACGCGTTGTTG-3'
 - b. Experimental primers:

pslA-F-qPCR: 5'-AAGATCAAGAAACGCGTGGAAT-3'

pslA-R-qPCR: 5'-TG TAGAGGT CGAACCACACCG-3'

pelA-F-qPCR: 5'-CCTTCAGCCATCCGTTCTTCT-3'

pelA-R-qPCR: 5'-TCGCGTACGAAGTCGACCTT-3'
4. Record information regarding your qPCR run according to the Minimal Information for Publication of Quantitative Real-Time PCR Experiments (MIQE) Guidelines (Bustin *et al.*, 2009 and 2010).

Data analysis

1. Experiments should be performed in biological triplicates. For this protocol, a biological replicate is defined as a biofilm grown, harvested, and flow-sorted on a separate experimental day.
2. We used the comparative Ct method ($\Delta\Delta C_t$) (Schmittgen and Livak, 2008) to determine the fold change in target gene expression between the reporter "on" and "off" populations.
 - a. From the qRT-PCR dataset, calculate the ΔC_t for each experimental gene (*e.g.*, *pelA*, *pslA*) in either the sorted "on" or "off" populations by subtracting the Ct of the control gene (*ampR*) from the experimental gene's Ct. For example, this is the ΔC_t for *pelA* in biofilm cells that were sorted as having high c-di-GMP reporter activity:

$$\Delta C_t = C_{t_{pelA-ON}} - C_{t_{ampR-ON}}$$

- b. Then, calculate the $\Delta\Delta C_t$ by subtracting the ΔC_t of *pelA* in the reporter "off" population from the ΔC_t of *pelA* in the reporter "on" population.
 - c. Finally, to determine the fold change in expression in the reporter "on" population compared to the reporter "off" population, take $2^{-\Delta\Delta C_t}$. Data can now be presented as the average fold change in *PelA* or *PslA* expression in the *P_{cdrA}* reporter sorted "on" population (high GFP)

relative to the “off” population (low GFP) for the three biological replicates (see Figure 4 for an example of how these data can be presented).




Figure 4. Average fold change in *pelA* and *pslA* mRNA transcript levels between the sorted c-di-GMP “off” and “on” populations. Levels of expression of Pel or Psl mRNA in the “on” population were normalized to the “off” population. This figure was reprinted from Figure 1E in Armbruster *et al.* (2019).

Recipes

1. LB Broth (Miller’s)

10 g Peptone

5 g Gibco Bacto Yeast extract

10 g NaCl

Ingredients should be combined in one liter of purified and deionized water, then autoclaved to sterilize.

2. Modified Jensen’s Chemically Defined Media (Jensen *et al.*, 1980; Jennings *et al.*, 2015)

*NaCl, 85.6 mM

*K₂HPO₄, 14.4 mM

*Sodium glutamate, 92 mM

*Valine, 24 mM

*Phenylalanine, 8 mM

Glucose, 70 mM

MgSO₄, 1.33 mM

CaCl₂, 0.14 mM

FeSO₄, 0.0039 mM

ZnSO₄, 0.0085 mM

Ingredients with an asterisk should be combined, pH adjusted to 7.3, and then autoclaved to sterilize. The remaining ingredients (metals and glucose) should be prepared separately and filter-sterilized. After autoclaving, allow the solution to cool to room temperature, then add the remaining filter-sterilized ingredients.

3. Phosphate-buffered saline (PBS), pH 7.4

8 g NaCl

0.2 g KCl

1.44 g Na₂HPO₄

0.24 g KH₂PO₄

Combine the ingredients to a final volume of 950 ml in purified and deionized water, adjust pH to 7.4, then add more water until the final volume is 1 L.

Acknowledgments

This work was supported by a Cystic Fibrosis Foundation Carol Basbaum Memorial Research Fellowship (ARMBRU19F0) and National Institutes of Health T32HL129949 to CRA and R01AI077628 and R01AI143916-01 to MRP. This protocol was used in a recently published work by the authors (Armbruster *et al.*, 2019). Illustrations were created with biorender.com.

Competing interests

The authors have no competing interests to declare.

Ethics

No human or animal specimens used.

References

1. Ackermann, M. (2015). [A functional perspective on phenotypic heterogeneity in microorganisms](#). *Nat Rev Microbiol* 13(8): 497-508.
2. Andersen, J. B., Sternberg, C., Poulsen, L. K., Bjorn, S. P., Givskov, M. and Molin, S. (1998). [New unstable variants of green fluorescent protein for studies of transient gene expression in bacteria](#). *Appl Environ Microbiol* 64(6): 2240-2246.
3. Armbruster, C. R., Lee, C. K., Parker-Gilham, J., de Anda, J., Xia, A., Zhao, K., Murakami, K., Tseng, B. S., Hoffman, L. R., Jin, F., Harwood, C. S., Wong, G. C. and Parsek, M. R. (2019). [Heterogeneity in surface sensing suggests a division of labor in *Pseudomonas aeruginosa* populations](#). *Elife* 8: e45084.
4. Arora, S. K., Ritchings, B. W., Almira, E. C., Lory, S. and Ramphal, R. (1997). [A transcriptional activator, FleQ, regulates mucin adhesion and flagellar gene expression in *Pseudomonas aeruginosa* in a cascade manner](#). *J Bacteriol* 179(17): 5574-5581.
5. Baraquet, C. and Harwood, C. S. (2016). [FleQ DNA Binding Consensus Sequence Revealed by Studies of FleQ-Dependent Regulation of Biofilm Gene Expression in *Pseudomonas aeruginosa*](#). *J Bacteriol* 198(1): 178-186.

6. Baraquet, C., Murakami, K., Parsek, M. R. and Harwood, C. S. (2012). [The FleQ protein from *Pseudomonas aeruginosa* functions as both a repressor and an activator to control gene expression from the *pel* operon promoter in response to c-di-GMP.](#) *Nucleic Acids Res* 40(15): 7207-7218.
7. Borlee, B. R., Goldman, A. D., Murakami, K., Samudrala, R., Wozniak, D. J. and Parsek, M. R. (2010). [Pseudomonas aeruginosa uses a cyclic-di-GMP-regulated adhesin to reinforce the biofilm extracellular matrix.](#) *Mol Microbiol* 75(4): 827-842.
8. Bustin, S. A., Benes, V., Garson, J. A., Hellemans, J., Huggett, J., Kubista, M., Mueller, R., Nolan, T., Pfaffl, M. W., Shipley, G. L., Vandesompele, J. and Wittwer, C. T. (2009). [The MIQE guidelines: minimum information for publication of quantitative real-time PCR experiments.](#) *Clin Chem* 55(4): 611-622.
9. Bustin, S.A., Beaulieu, J. F., Huggett, J., Jaggi, R., Kibenge, F.S., Olsvik, P.A., Penning, L.C. and Toegel, S. (2010). [MIQE précis: Practical implementation of minimum standard guidelines for fluorescence-based quantitative real-time PCR experiments.](#) *BMC Mol Biol* 11: 74.
10. Güvener, Z. T. and Harwood, C. S. (2007). [Subcellular location characteristics of the *Pseudomonas aeruginosa* GGDEF protein, WspR, indicate that it produces cyclic-di-GMP in response to growth on surfaces.](#) *Mol Microbiol* 66(6): 1459-1473.
11. Hengge, R. (2009). [Principles of c-di-GMP signalling in bacteria.](#) *Nat Rev Microbiol* 7(4): 263-273.
12. Hickman, J. W. and Harwood, C. S. (2008). [Identification of FleQ from *Pseudomonas aeruginosa* as a c-di-GMP-responsive transcription factor.](#) *Mol Microbiol* 69(2): 376-389.
13. Hickman, J. W., Tifrea, D. F. and Harwood, C. S. (2005). [A chemosensory system that regulates biofilm formation through modulation of cyclic diguanylate levels.](#) *Proc Natl Acad Sci U S A* 102(40): 14422-14427.
14. Huangyutham, V., Guvener, Z. T. and Harwood, C. S. (2013). [Subcellular clustering of the phosphorylated WspR response regulator protein stimulates its diguanylate cyclase activity.](#) *mBio* 4(3): e00242-00213.
15. Jennings, L.K., Storek, K.M., Ledvina, H.E., Coulon, C., Marmont, L.S., Sadovskaya, I., Secor, P.R., Tseng, B.S., Scian, M., Filloux, A., et al. (2015). [Pel is a cationic exopolysaccharide that cross-links extracellular DNA in the *Pseudomonas aeruginosa* biofilm matrix.](#) *Proc Natl Acad Sci U S A* 112: 11353-11358.
16. Jensen, S. E., Fecycz, I. T. and Campbell, J. N. (1980). [Nutritional factors controlling exocellular protease production by *Pseudomonas aeruginosa*.](#) *J Bacteriol* 144(2): 844-847.
17. Kulasakara, H., Lee, V., Brencic, A., Liberati, N., Urbach, J., Miyata, S., Lee, D. G., Neely, A. N., Hyodo, M., Hayakawa, Y., Ausubel, F. M. and Lory, S. (2006). [Analysis of *Pseudomonas aeruginosa* diguanylate cyclases and phosphodiesterases reveals a role for bis-\(3'-5'\)-cyclic-GMP in virulence.](#) *Proc Natl Acad Sci U S A* 103(8): 2839-2844.

18. Rybtke, M. T., Borlee, B. R., Murakami, K., Irie, Y., Hentzer, M., Nielsen, T. E., Givskov, M., Parsek, M. R. and Tolker-Nielsen, T. (2012). [Fluorescence-based reporter for gauging cyclic di-GMP levels in *Pseudomonas aeruginosa*](#). *Appl Environ Microbiol* 78(15): 5060-5069.
19. Schmittgen, T. D. and Livak, K. J. (2008). [Analyzing real-time PCR data by the comparative C\(T\) method](#). *Nat Protoc* 3(6): 1101-1108.
20. Zaitoun, I., Erickson, C.S., Schell, K. and Epstein, M. L. (2010). [Use of RNAlater in fluorescence-activated cell sorting \(FACS\) reduces the fluorescence from GFP but not from DsRed](#). *BMC Res Notes* 3: 328-328.