

Preparation and Characterization of Lipid Nanoparticles Co-loaded With DNA and Nitro-Oleic Acid

Manthan N. Patel^{1,2,*}, Sachchidanand Tiwari^{2,3}, and Jacob S. Brenner^{1,2,*}

¹Department of Systems Pharmacology and Translational Therapeutics, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA, USA

²Division of Pulmonary, Allergy, and Critical Care, Department of Medicine, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA, USA

³Department of Microbiology, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA, USA

*For correspondence: pmanthan@gmail.com; jacob.brenner@pennmedicine.upenn.edu

Abstract

Lipid nanoparticles (LNPs) are powerful carriers for nucleic acid delivery, but plasmid DNA-loaded LNPs (pDNA-LNPs) have been limited by inflammation and toxicity. We showed that standard pDNA-LNPs activate the cGAS–STING pathway, leading to severe immune responses and mortality in mice. To overcome this, we co-loaded nitro-oleic acid (NOA), an endogenous STING inhibitor, into pDNA-LNPs. NOA-pDNA-LNPs mitigated inflammation, enabled safe in vivo delivery, and supported sustained gene expression for months. Here, we present a detailed protocol for producing and characterizing NOA-pDNA-LNPs to facilitate safer, long-term gene delivery applications.


Key features

- Provides a step-by-step protocol to produce plasmid DNA–LNPs co-loaded with nitro-oleic acid (NOA), optimized for both in vitro and in vivo applications.
- Includes methods for quantitative assessment of DNA and NOA encapsulation efficiencies, particle size, and quality control metrics.

Keywords: Lipid nanoparticles, Non-viral, Drug delivery, pDNA, Anti-inflammatory lipids, cGAS-STING

This protocol is used in: Nat Biotechnol (2025), DOI: 10.1038/s41587-025-02556-5

Graphical overview



Lipid nanoparticle (LNP) production process. 1. Prepare lipid and DNA mixes in separate 2 mL Eppendorf tubes. 2. Mix the two phases using either a microfluidic chip or a vortex mixer. 3. Dialyze the formed LNPs against 1× PBS, pH 7.4, for 2 h to remove remaining ethanol and adjust pH. 4. Validate LNP quality before use or storage.

Background

The success of COVID-19 mRNA vaccines showcased the power of lipid nanoparticles (LNPs) to deliver nucleic acids and achieve high protein expression [1]. While most LNP applications focus on mRNA or siRNA, these cargos are short-lived and lack promoter control. In contrast, DNA offers longer-term expression (months to years), cell type-specific promoter regulation, and greater stability, making it ideal for expressing therapeutic proteins, gene editors, or shRNAs for long-term effects [2–4].

DNA-loaded LNPs combine durable expression with the advantages of LNPs: high delivery efficiency, low immunogenicity (compared to viral vectors), and flexible cargo capacity. However, their development has lagged due to safety concerns. We found that plasmid DNA-LNPs (pDNA-LNPs) trigger strong inflammation in naive mice, leading to 100% mortality at standard doses [5]. Using knockout mouse models and *in vitro* assays, we identified the cGAS–STING pathway as the main driver. This innate immune pathway detects cytosolic DNA and activates pro-inflammatory cytokines like type I interferons and IL-6.

To overcome this, we took inspiration from DNA viruses that evade STING-mediated immunity. We loaded pDNA-LNPs with nitro-oleic acid (NOA), an endogenous STING inhibitor. These NOA-pDNA-LNPs eliminated inflammation *in vitro* and prevented mortality *in vivo*, without compromising transgene expression.

Together, these results demonstrate that pDNA-LNPs co-loaded with bioactive lipids like NOA can enable safe, long-term gene expression. Here, we present a detailed protocol for producing NOA-pDNA-LNPs to support a broader application of this promising gene delivery platform.

Materials and reagents

Reagents

1. SM-102 (Echelon Biosciences, catalog number: N-1020); store at -20 °C
2. Cholesterol (Avanti Polar Lipids, catalog number: 700100P); store at -20 °C
3. 1,2-distearoyl-sn-glycero-3-phosphocholine (DSPC) (Avanti Polar Lipids, catalog number: 850365); store at -20 °C
4. 1,2-dimyristoyl-rac-glycero-3-methoxypolyethylene glycol-2000 (DMG-PEG 2000) (Avanti Polar Lipids, catalog number: 880151); store at -20 °C

5. Nitro-oleic acid (NOA) (Echelon Biosciences, catalog number: L-0112); store at -20 °C
6. pALD-CV77-Luciferase (pDNA) (Aldevron, catalog number: 5078-5); store at 4 °C short-term and at -20 °C long-term
7. 200 proof ethanol (100% ethanol) (Millipore Sigma, catalog number: 459836-100ML); store at room temperature
8. DEPC water (Invitrogen, catalog number: AM9916); store at 4 °C
9. Phosphate-buffered saline (1× PBS), filter-sterilized (Thermo Fisher Scientific, catalog number: J61196.AP)
10. Acetonitrile (HPLC grade) (Fisher Scientific, catalog number: A998SK-1)
11. Trifluoroacetic acid (TFA) (Millipore Sigma, catalog number: 8082600026)
12. Quant-iT PicoGreen dsDNA assay (Invitrogen, catalog number: P11496); store at 4 °C short-term and at -20 °C long-term
13. Sodium citrate dihydrate (Millipore Sigma, catalog number: W302600)
14. Citric acid (Millipore Sigma, catalog number: 251275)
15. TE buffer (Thermo Fisher Scientific, catalog number: 12090015)

Laboratory supplies

1. 2 mL Eppendorf tubes (Fisher Scientific, catalog number: 05-402-7)
2. 15 mL conical Tubes (Thermo Fisher Scientific, catalog number: 339650)
3. 1–300 µL pipette tips (Fisher Scientific, SureOne™ Micropoint Pipette Tips, catalog number: 02-707-410)
4. 500 mL beaker (Corning, catalog number: 1000J-500)
5. Dialysis tubing 12–14 kD (Repligen, catalog number: 132697)
6. Dialysis tubing 55 mm and 35 mm closures (Repligen, catalog numbers: 132734 and 132760, respectively)
7. Disposable cuvettes (GMBH + Co KG, catalog number: 759075D)
8. Zeba columns (Thermo Fisher Scientific, catalog number: 89877)
9. UPLC vials (Thermo Fisher Scientific, catalog number: 6ESV9-04PP)

Equipment

1. P20 pipette (Gilson, catalog number: F144056M)
2. P200 pipette (Gilson, catalog number: F144058M)
3. Ultrasonic bath (Branson, catalog number: BRANCPX-952-538R)
4. NanoAssemblr™ Ignite™ instrument (Cytiva, catalog number: NIN0001)
5. NanoAssemblr™ Ignite™ cartridges (Cytiva, catalog number: NIN0061)
6. Vortex-Genie 2 (Fisher Scientific, catalog number: 50-728-002)
7. Zetasizer Pro ZS (Malvern Panalytical)
8. MiniMag Magnetic Stirrer (NETA/Sigma-Aldrich, catalog number: Z742293)
9. Acquity I-Class UPLC system (Waters)
10. Cortecs C18 column (1.6 µm particle size, 2.1 × 50 mm)

Procedure

A. Preparation of lipid mix

1. Thaw lipids to room temperature prior to opening the vial.
2. Dissolve the lipids using 100% ethanol to the indicated stock concentrations (Table 1).
 - a. Follow the pipetting protocol at <https://www.addgene.org/protocols/pipetting/>.
 - b. Some lipids may require brief heating to 50–60 °C and sonication for complete dissolution.

Caution: Ensure cholesterol crystals are completely dissolved.

- c. Store at -20 °C after aliquoting the stock solutions.
- d. This protocol is designed for a 40:1 w/w ratio of total lipids to DNA (excluding NOA in the total lipid weight).
 - i. 25 mM lipid mix for microfluidics.
 - ii. 6.25 mM lipid mix for vortex mixing.

3. Prepare a 2 mL Eppendorf tube and pipette the indicated volume of each lipid/ethanol (Table 1).

Notes:

1. Vortex mixing is typically used for producing LNPs for in vitro experiments, and microfluidics is used for in vivo experiments.
2. The “L” in NOA:L indicates the total lipid content excluding NOA.

Table 1. Components in the lipid mix tube

Component	Name	Stock concentration (mg/mL)	Molar (%) in LNP formulation	Volume (μL) for 0.4 mL of LNPs (microfluidics)	Volume (μL) for 0.16 mL of LNPs (vortex)
Ionizable lipid	SM-102	50	50	35.5	4.4
Cholesterol	Cholesterol	11.6	38.5	64.2	8.0
Helper lipid	DSPC	10	10	19.8	2.5
PEG lipid	DMG-PEG 2000	10	1.5	18.8	2.4
Ethanol	-	-	-	37.5	79.69
Anti-inflammatory lipid	NOA	13.5	0.2 NOA:L (mol:mol)	24.2	3.03

B. Preparation of DNA mix

1. Dissolve DNA to 1 mg/mL using nuclease-free water (TE buffer works as well).

2. Prepare a 2 mL Eppendorf tube and pipette the indicated volume of DNA and buffer (Table 2).

Note: To prepare a 50 mM citrate buffer at pH 4, dissolve 496 mg of sodium citrate dihydrate and 636 mg of citric acid in 100 mL of water. Adjust pH using HCl or NaOH if needed.

Table 2. Components in the DNA mix tube

Component	Name	Composition	Volume (μL) for 0.4 mL of LNPs (microfluidics, in vivo)	Volume (μL) for 0.16 mL of LNPs (vortex, in vitro)
DNA	pALD-CV77-Luciferase	1 mg/mL in water	50.7	3.9
Buffer	Citrate buffer	50 mM, pH 4	348.3	116.1

C. Production of LNPs using microfluidics

1. Insert the microfluidic cartridge into the NanoAssemblr™ Ignite™ (Figure 1).



Figure 1. Microfluidic chip for NanoAssemblr™ Ignite™

2. Load one 1 mL BD syringe with the lipid mix and another with the DNA mix (Figure 2).

Tip: To remove bubbles, invert the syringe and gently tap/flick until the bubble flows to the top. Then, push the liquid until it reaches the top of the tip.

Critical: Ensure the volumes are > 0.3 mL in the DNA syringe and > 0.1 mL in the lipid syringe.

3. Add waste (left) and LNP collection (right) 15 mL conical tubes to the right side of the cartridge/syringes (Figure 2).




Figure 2. Use 1 mL BD syringes to load lipid and DNA mixes into the microfluidic chip

4. Enter various formulation parameters (Figure 3):

- a. Syringe type: BD 1 mL.
- b. Flow rate ratio: 3:1 (DNA to lipid).
- c. Total LNP volume: 0.4 mL.
- d. Total flow rate: 12 mL/min.
- e. Waste: 0.02 mL for start and end waste.




Figure 3. Lipid nanoparticle (LNP) formulation parameters

Note: A detailed protocol and machine function/troubleshooting can be found in the NanoAssemblr™ Ignite™ manual.

5. Continue to Section E after the machine finishes its run.

D. Production of LNPs using vortex mixing

1. Pipette the various lipids into a 2 mL Eppendorf tube (Table 1).

Tip: Add ethanol first as it has the highest volume.

2. Pipette the DNA and citrate buffer into another 2 mL Eppendorf tube (Table 2).

3. While vortexing the DNA mix tube at maximum speed (setting 10 on the Vortex-Genie 2), quickly pipette 40 μ L of the lipid mix all at once (Figure 4). The pipette tip should be positioned just above the surface of the DNA mix, not submerged.

4. Continue to vortex for 10 s after adding the lipid mix, then proceed to Section E.



Figure 4. Vortex mixing for lipid nanoparticle (LNP) production

E. Dialysis of LNPs (Figure 5)

1. Cut and pre-wet the dialysis tubing (10 kDa MW cutoff) for 1 min.

2. Clip the bottom end of the tubing using the 35 mm weighted clip and pipette in LNPs.

Tip: Squeeze the tubing from the sides first for easier LNP loading.

3. Clip the top end of the tubing using a 55 mm unweighted clip to secure the LNPs inside the tubing.

4. Submerge in 1 \times PBS, pH 7.4, and stir for 2 h at room temperature.

Caution: Ensure the LNP solution is completely submerged in the dialysis bath.




Figure 5. Dialyzing lipid nanoparticles (LNPs) against 1 \times PBS. (A) Pre-wetted dialysis tubing. (B) LNPs loaded in the

tubing with both ends clipped. (C) Dialysis of LNPs with stirring.

5. After dialysis, remove the top clip and pipette out LNPs into a fresh 2 mL Eppendorf tube.

Pause point: Store the LNPs at 4 °C or proceed to Section F.

Caution: It is best practice to use the LNPs within a couple of days.

F. Characterization of LNPs

F1. Measure particle size using dynamic light scattering (DLS)

1. Dilute the LNPs in 1× PBS solution directly in the disposable DLS cuvettes.

a. For vortex-made LNPs: 20 µL of LNPs into 480 µL of 1× PBS.

b. For microfluidic-made LNPs: 5 µL of LNPs into 495 µL of 1× PBS.

2. Insert in the Zetasizer Pro ZS and conduct size and PDI measurement of the LNPs.

Note: A detailed protocol can be found in the Zetasizer Pro ZS (Malvern) manual.

F2. Measure DNA encapsulation

1. Follow the exact protocol from the Quant-iT™ PicoGreen™ dsDNA Assay kits and dsDNA reagents from ThermoFisher Scientific.

2. Typical LNP dilutions:

a. For vortex LNPs, dilute 40 µL of LNPs in 210 µL of TE buffer.

b. For microfluidic LNPs, dilute 10 µL of LNPs in 240 µL of TE buffer.

F3. Measure NOA encapsulation

1. Purify out unencapsulated NOA using Zeba size exclusion columns.

Note: A detailed protocol can be found in the Zeba column manual (Thermo Fisher).

a. Remove Zeba column storage solution by centrifuging at 1,500× g for 1 min.

b. Wash the column using 300 µL of 1× PBS, three times, centrifuging at 1,500× g for 1 min between each wash.

c. Load 70 µL of LNPs and collect in a fresh 2 mL Eppendorf tube.

2. Dilute LNPs (before Zeba column purification and after purification) 10× into 2 mL polypropylene screw-top microvials containing 50% water and 50% acetonitrile: 10 µL of LNP into 45 µL of water + 45 µL of acetonitrile. Briefly vortex to mix.

3. Run the two samples on an Acquity I-Class UPLC system with Cortecs C18 column (1.6 µm particle size, 2.1 × 50 mm).

a. Run with mobile phase A [water with 0.1 vol% trifluoroacetic acid (TFA)] and mobile phase B (acetonitrile with 0.1 vol% TFA).

Caution: TFA is a strong, corrosive, and volatile acid

b. Set the flow rate to 0.5 mL/min with a gradient of 10% B from 0 min to 0.4 min, linearly increasing to 90% B from 0.4 min to 0.8 min, holding at 90% B from 0.8 min to 1.3 min, and then linearly decreasing to 10% B from 1.3 min to 1.4 min.

c. Set the sample injection volume to 8 µL.

4. Using a standard curve for NOA, calculate the total concentration of NOA in LNP.

5. Based on the pre-Zeba and post-Zeba NOA concentrations, determine the NOA encapsulation efficiency.

Data analysis

1. DLS analysis showed that the average size for NOA-pDNA-LNPs (Figure 6) is 80.7 ± 2.0 nm, based on triplicate measurements.

Note: LNPs made via the vortex method typically are larger than those made with microfluidics. Inefficient mixing may lead to aggregation of LNPs (Figure 7). The DNA encapsulation efficiency and transfection efficiency may be slightly lower as well.




Figure 6. Size measurement of nitro-oleic acid (NOA)-plasmid DNA (pDNA)-loaded lipid nanoparticles (LNPs) using dynamic light scattering (DLS)




Figure 7. An example of a rejected lipid nanoparticle (LNP) sample. Arrows indicate the presence of aggregates.

2. The PicoGreen analysis showed that the average pDNA encapsulation for NOA-pDNA-LNPs is $87.7\% \pm 1.67\%$, based on triplicate measurements.
3. The ultra-performance liquid chromatography (UPLC) analysis showed that the average NOA encapsulation for NOA-pDNA-LNPs is $82.8\% \pm 1.72\%$, based on triplicate measurements.

Note: Elution of NOA (1.5 min, wavelength: 262 nm) is depicted in Figure 8 when using the gradient method, as outlined in Section F3. This elution time should be similar for both vortex and microfluidic LNPs.




Figure 8. Ultra-performance liquid chromatography (UPLC) chromatogram of nitro-oleic acid (NOA) content within NOA-pDNA-lipid nanoparticles (LNPs)

Validation of protocol

This protocol (or parts of it) has been used and validated in the following research article:

- Patel et al. [5]. Safer non-viral DNA delivery using lipid nanoparticles loaded with endogenous anti-inflammatory lipids. Nature Biotechnology (Figure 3b).

General notes and troubleshooting

General notes

1. Batch-to-batch variability in LNP properties (size, encapsulation efficiency) may occur due to slight differences in mixing technique, lipid lot quality, or ambient temperature. It is recommended to validate each new lipid lot and calibrate the mixing process when switching operators.
2. This protocol was optimized for the delivery of pDNA using the pALD-CV77-Luciferase plasmid (Aldevron). While it may be applicable to other plasmids or DNA types (e.g., minicircle, dbDNA, or ssDNA), optimization of lipid ratios and buffer conditions may be required.
3. The stability of LNPs is limited. It is best to use freshly prepared LNPs within 2–3 days of production. Extended storage at 4 °C may reduce transfection efficiency due to LNP aggregation or hydrolysis.
4. pDNA-LNPs containing NOA show significantly reduced inflammation compared to standard LNPs. However, further validation is recommended before clinical or immunologically sensitive applications.
5. Downstream quantification assays (e.g., PicoGreen and UPLC for NOA) are sensitive to contamination. Use nuclease-free and LC/MS-grade reagents, clean pipettes, and avoid repeated freeze-thaw cycles of standards or reagents.
6. NanoAssemblr™ Ignite™ provides reproducible LNPs, but settings may require tuning when switching to different nucleic acids, buffers, or lipid compositions. Refer to the manufacturer's manual and consider pilot runs for each new formulation.
7. Both vortex- and microfluidic-produced LNPs can be used for in vitro and in vivo studies, though microfluidic LNPs are preferred for better reproducibility.
8. Typical acceptance criteria are as follows: size < 200 nm; PDI < 0.25; and DNA encapsulation efficiency > 80%. Failure in one or more of these criteria warrants rejection of the LNP batch.
9. This protocol is amenable to different pDNA sizes (tested for 3–7 kb) and different NOA percentages (tested for 0.1–1 NOA/L).

Troubleshooting

Problem 1: Low LNP yield after microfluidic mixing.

Possible cause: Incomplete loading of syringes or improper flow rate settings.

Solution: Ensure DNA and lipid syringe volumes exceed 0.3 and 0.1 mL, respectively, and verify if flow rate settings on the NanoAssemblr™ Ignite are correctly entered (TFR: 12 mL/min, FRR: 3:1).

Problem 2: Cloudy or phase-separated LNP solution.

Possible cause: Incomplete lipid dissolution before mixing.

Solution: Briefly heat and sonicate lipid stock (especially cholesterol) to ensure full solubilization before preparing the ethanol lipid mix.

Problem 3: Large or polydisperse particle sizes in DLS analysis.

Possible cause: Improper mixing technique or ethanol not fully removed post-dialysis.

Solution: Verify if a correct and vigorous vortexing or microfluidic mixing was used. Ensure a full 2-h dialysis in 1× PBS and that the LNPs are fully submerged and stirred.

Problem 4: Low DNA encapsulation efficiency.

Possible cause: Incorrect lipid-to-DNA ratio or degraded DNA.

Solution: Confirm 40:1 w/w lipid:DNA ratio was used, and that the DNA stock is intact (check via gel electrophoresis if needed). Use freshly prepared or properly stored DNA at 1 mg/mL.

Problem 5: Difficulty loading LNPs into dialysis tubing.

Possible cause: Tubing was not pre-wet, or tubing collapsed during loading.

Solution: Pre-wet tubing with water before use, and pinch sides gently to open the channel for easier pipetting.

Problem 6: No peak or inconsistent NOA signal on UPLC.

Possible cause: Incorrect gradient or injection parameters; TFA evaporation or contamination.

Solution: Double-check mobile phase composition and gradient timing. Prepare fresh mobile phases with 0.1% TFA.

Problem 7: Bubbles in syringes during microfluidic mixing.

Possible cause: Air trapped during syringe filling.

Solution: Invert syringes and gently flick to move bubbles to the top. Push the plunger slowly to expel bubbles before insertion into the NanoAssemblr™.

Acknowledgments

Research reported in this publication was supported by the American Heart Association under grant 24PRE1195406 (to M.N.P.), grants R01-HL-153510, R01-HL160694, and R01-HL164594 (to J.S.B.), and grant R01AI153064 (to N.P.). This protocol was used in [5].

Competing interests

M.N.P., S.T., and J.S.B. have filed a patent for the DNA-LNP technology described herein. The other authors declare no competing interests.

Received: July 15, 2025; Accepted: August 20, 2025; Available online: September 01, 2025; Published: September 20, 2025

References

1. Hou, X., Zaks, T., Langer, R. and Dong, Y. (2021). Lipid nanoparticles for mRNA delivery. *Nat Rev Mater.* 6(12): 1078–1094. <https://doi.org/10.1038/s41578-021-00358-0>
2. Hughes, T. S., Langer, S. J., Johnson, K. W., Chavez, R. A., Watkins, L. R., Milligan, E. D. and Leinwand, L. A. (2009). Intrathecal Injection of Naked Plasmid DNA Provides Long-term Expression of Secreted Proteins. *Mol Ther.* 17(1): 88–94. <https://doi.org/10.1038/mt.2008.230>
3. Zhu, Y., Shen, R., Vuong, I., Reynolds, R. A., Shears, M. J., Yao, Z. C., Hu, Y., Cho, W. J., Kong, J., Reddy, S. K., et al. (2022). Multi-step screening of DNA/lipid nanoparticles and co-delivery with siRNA to enhance and prolong gene expression. *Nat Commun.* 13(1): e1038/s41467-022-31993-y. <https://doi.org/10.1038/s41467-022-31993-y>
4. Chen, Z. Y., He, C. Y., Ehrhardt, A. and Kay, M. A. (2003). Minicircle DNA vectors devoid of bacterial DNA result in persistent and high-level transgene expression in vivo. *Mol Ther.* 8(3): 495–500. [https://doi.org/10.1016/s1525-0016\(03\)00168-0](https://doi.org/10.1016/s1525-0016(03)00168-0)
5. Patel, M. N., Tiwari, S., Wang, Y., O'Neill, S., Wu, J., Omo-Lamai, S., Espy, C., Chase, L. S., Majumder, A., Hoffman, E., et al. (2025). Safer non-viral DNA delivery using lipid nanoparticles loaded with endogenous anti-inflammatory lipids. *Nat Biotechnol.*: e1038/s41587-025-02556-5. <https://doi.org/10.1038/s41587-025-02556-5>