

## Bioassay of Extracts of the Endophytic Fungi

Susheel Kumar<sup>1</sup> and Nutan Kaushik<sup>2\*</sup>

<sup>1</sup>Bioresources and Biotechnology, The Energy and Resources Institute, New Delhi, India;

<sup>2</sup>Environmental and Industrial Biotechnology, The Energy and Resources Institute, New Delhi, India

\*For correspondence: kaushikn@teri.res.in

[Abstract] Many of the microbes including fungi produce metabolites possessing antifungal activity and in some cases fungal metabolites are main cause of antifungal activity resulted by fungus. To infer that antifungal activity is due to fungal metabolites, there is a need to develop a repeatable procedure to assay these metabolites. Here we are presenting the poisoned food technique for bioassay of extract from endophytic fungi, where culture media is supplemented with extract to testpathogen viability as proxy for antifungal activity.

## **Materials and Reagents**

- 1. Culture of plant pathogenic fungi (Sclerotinia sclerotiorum)
- 2. Methanol (Qualigens)
- 3. Potato dextrose agar (PDA) media (HiMedia Laboratories)
- 4. Commercial disinfectant( 70%ethanol)
- 5. Distilled water
- 6. PDA media (see Recipes)

## **Equipment**

- 1. Measuring cylinder
- 2. Inoculating needle
- 3. Spirit lamp
- 4. Cork borer
- 5. Conical flask
- 6. Measuring scale
- 7. 10 cm Petri dishes
- 8. pH meter (Eutech Instruments pH tutor)
- 9. Autoclave (Nat steel)
- 10. BOD incubator (Toshiba)



- 11. Laminar air flow hood (Toshiba)
- 12. Micropipette (Eppendorf)
- 13. Microbalances (Sartorious, model: RC210P)

# **Procedure**

- 1. Solvent extracts from endophytic fungi were tested against *Sclerotinia sclerotiorum* a well known plant pathogenic fungus.
- 2. 30 mg of dried extract was dissolved in 800 µl of methanol.
- 3. 200 and 400  $\mu$ l of the dissolved extract were added to 30 ml of molten PDA media, mixed and then poured equally into three 10 cm Petri plates to obtain 250 and 500  $\mu$ g/ml extract concentrations, respectively.
- 4. To obtain 1,000  $\mu$ g/ml concentration, 30 mg of dried extracts were dissolved in 400  $\mu$ l of methanol and 400  $\mu$ l were added to 30 ml of molten PDA.
- 5. Control growth plates contained 400 µl of methanol.
- 6. *S. sclerotiorum* was inoculated at the centre of the plate and radial growth was measured at 48 h intervals till the control plate attained the full growth.
- 7. Per cent growth inhibitions (%GI) of the extracts were calculated relative to the growth on the control plate.

 $%GI = \{(A-B)/A\} \times 100$ 

Where A = Growth of plant pathogenic fungus in control plate

B = Growth of plant pathogenic fungus plates supplemented with extract from endophytic fungus

#### Recipes

1. PDA media

39.1 g of PDA media was dissolved in 1,000 ml of distilled water and pH was adjusted to 7.0

## **Acknowledgments**

This protocol was adopted from Kumar and Kaushik (2013). Authors are grateful to their host institution, The Energy and Resources Institute (TERI), New Delhi, India for funding the research. Susheel Kumar is grateful to University Grant Commission, New Delhi for a research fellowship.



# References

1. Kumar, S. and Kaushik, N. (2013). <u>Endophytic fungi isolated from oil-seed crop Jatropha curcas produces oil and exhibit antifungal activity</u>. *PLoS One* 8(2): e56202.